ULTRASONIC CLEANING MONITOR

For Use in Routine Monitoring of the Wash Phase
in Ultrasonic Cleaners.

INDICATIONS FOR USE:

The single use Hu-Friedy Ultrasonic Cleaning Monitor is designed for monitoring cleaning efficiency of ultrasonic cleaners.

RECOMMENDED USE:

Daily monitoring of ultrasonic cleaners for machine release, and periodically throughout the day for load release.

DAILY MACHINE RELEASE

1. Insert the Ultrasonic Cleaning Monitor strip (IMS-1200U) into the Cleaning Monitor Holder (IMS-1200H). Insert so that the BLUE SQUARE is fully exposed.
2. Place the Ultrasonic Cleaning Monitor and holder in the ultrasonic cleaner with cleaning liquid but no instruments. Make sure the bluesquare faces up.
3. Run a normal ultrasonic cleaner cycle.
4. At the end of the cycle, remove the Ultrasonic Cleaning Monitor and holder from the ultrasonic cleaner and review the results. If the cleaning cycle was effective, the BLUE SQUARE of test soil should disappear.
5. If desired, staple the Ultrasonic Cleaning Monitor strip onto the Ultrasonic Cleaning Monitor Record Sheet and complete the appropriate information. Retain as a permanent record.

ROUTINE LOAD RELEASE

1. Insert the Ultrasonic Cleaning Monitor strip into the Cleaning Monitor Holder. Insert so that the BLUE SQUARE is fully exposed.
2. Place the Ultrasonic Cleaning Monitor and holder in the ultrasonic cleaner on top of the instruments. Make sure the blue square is face up.
3. Process the ultrasonic cleaning load per normal procedure.
4. At the end of the cycle, remove the Ultrasonic Cleaning Monitor and holder from the ultrasonic cleaner and review the results. If the cleaning cycle was effective, the BLUE SQUARE of test soil should disappear.
5. If desired, staple the Ultrasonic Cleaning Monitor strip onto the Ultrasonic Cleaning Monitor Record Sheet and complete the appropriate information. Retain as a permanent record.

INTERPRETATION OF RESULTS

Examine the test strip for presence of any blue test soil remaining. In an efficient ultrasonic cycle, the blue square area should disappear. The presence of a large amount of blue test soil remaining on the device means a serious deficiency in ultrasonic cleaner efficiency. A small amount of blue remaining indicates that conditions are slightly deficient (marginal).

Untreated Monitor, or critical/poor result
(entire blue square is visible)

Insufficient result (some blue remaining in square)

Optimum result (no blue remaining in square)

STORAGE AND HANDLING

- Store Cleaning Monitors inside the protective foil bag provided at 40-86°F (4-30°C). Do not remove the monitors from the foil bag until ready to use.
- Keep away from disinfectants and fumes.

EXPIRY DATE

The expiry date is printed on the product packaging.

Please see website for Ultrasonic Cleaning Monitor Record Sheet.

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