



# **HU-FRIEDY** STERILISATION CONTAINER SYSTEM

Hu-Friedy's Sterilisation Containers are a reusable, rigid system used for the packaging of cassettes and instruments during sterilisation. An alternative to the wrap and pouch materials commonly used to package cassettes and instruments, Sterilisation Containers provide significant economic and environmental benefits along with enhanced safety.



- Two choices of filters to meet customers' needs
- Enhanced locking device to secure the container
- Compatible with all B-autoclaves
- · Improved filter support for a secure fit
- Lid and bottom are made of anodized aluminum making them lightweight with excellent drying properties. The aluminum alloy with an anodized oxide surface helps resist corrosion
- Two external display windows provide space for a process indicator card and container identification



### Please visit us at www.Hu-Friedy.eu



INSTRUMENT MANAGEMENT

# **HU-FRIEDY** STERILISATION CONTAINER SYSTEM

# A valid alternative to the traditional wrapping materials.



#### FOR 1 CASSETTE

# **IMDINCOOM**

**DIN Container for** max. 1 Container Cassette Top and bottom with filter\*

Dimensions: 312 x 190 x 50 mm.

\* Permanent (PTFE) or Disposable Paper filters to be ordered separately.



# FOR 1 CASSETTE

## **IMDINCO1M DIN Container for** max. 1 Container Cassette Top and bottom with filter\* Dimensions: 312 x 190 x 65 mm.



#### FOR 2 CASSETTES

IMDINCO2M **DIN Container for** max. 2 Container Cassettes Top and bottom with filter\* Dimensions: 312 x 190 x 100 mm.



#### FOR 3 CASSETTES

**IMDINCO3M DIN Container for** max. 3 Container Cassettes Top and bottom with filter\* Dimensions: 312 x 190 x 135 mm.

# IMS INFINITY CONTAINER CASSETTE

Fits for all DIN sterilisation containers. With handle on lid and filter outlets.

FOR 16 INSTRUMENTS

# **IMECC12 S**

Surgical container cassette for 16 Instruments, 2 clips for hinged instruments and accessory space

**Dimensions:** 268 x 176 x 34 mm



# **FOR 12 INSTRUMENTS** IMECC12

For 12 struments, with an accessory space Dimensions: 268 x 176 x 34 mm.



#### **CASSETTE COLOR CODING**

- 0 = Red 1 = Grev
- 2 = White
- 3 = Orange
- 4 = Purple
- 5 = Yellow 6 = Lavender
- 7 = Black
- 8 = Blue
- 9 = Green
- 11 = Ocean Blue

# **CONTAINER ACCESSORIES**

# (Artiliado) 223 ČE R

(11)

Statistics of the local division of the

And Berner, 1994, 50

\* without CE

•

#### **IMCO-FIL1M** Permanent

filters (PTFE) 2 units. Usable for up to 1200 sterilisation cvcles **Dimensions:** 



# **IMCO-FIL2M**



100 units per package

**Dimensions:** 118 x 235 mm



## **IMCO-IND** Indicators/Labels 100 units.



**IMCOM-SEAL** Sealing\* 100 units. Single use



### IMCO-C

Identification plates\* If with engraving, please indicate desired engraving when ordering

Color Codes : Red - IMCO-CO Orange - IMCO-C3 Gold - IMCO-C4 Blue - IMCO-C8 Green - IMCO-C9



# **IMS-COBL**

For reading sterilisation containers and sterilisation inventory. HIBC UIM-Standard, label with guide rail, attachment in standard slot Material: PPSU white, carbonizable. Designed to withstand type B steam sterilisation up to 140°C. Neutral (without Text)

#### Dimensions: 50 x 18 mm., 4 mm. thick

The barcode label is also available for the cassette Partcode: IMS-CABL

Manufacturer: Hu-Friedy Mfg. Co., LLC | 3232 N. Rockwell Street | Chicago, IL 60618 | USA Customer Care Department: Hu-Friedy Mfg. Co., LLC. Free Call +800 HUFRIEDY (00800 48 37 43 39) | Free Fax 00800 48 37 43 40 E-Mail: info@hu-friedy.eu | Website: Hu-Friedy.eu

©2017 Hu-Friedy Mfg. Co., LLC. All rights reserved. The contents of this Product Catalog & Reference Guide including, but not limited to, all text, photographs and graphics, are the sole and exclusive property of Hu-Friedy Mfg. Co., LLC and may not be copied, reproduced or distributed without prior written permission. HF-GB-CON/0317 Hu-Friedy Mfn Co. LLC reserves the right to change features of products images and text at any time without notice

How the best perform

