



IMS STERILISATION WRAP SOLUTIONS

IMS MONITOR TAPE

IMS Bonded Double Wrap

- Two sheets of wrap ultrasonically bonded together on two edges for quick and easy wrapping application, assuring compliance.
- Two different colors sheets to easily identify breaches post sterilization.
- High strength and durability aid in resistance to tears and punctures, protecting contents of the package during storage and transportation.
- Made of SMS (Spunbond Meltblown Spunbond) polypropylene material – any unused or unsaturated wrap can be recycled as a recycle code of "5".
- Compatible with many common sterilization cycles.
- Product part code and size printed on each sheet for quick and easy identification when reordering product.
- Not made with natural rubber latex, enhances safety for patients and users who suffer latex allergies.



IMS Sterilization Wrap

IMS Cassette Wrap is engineered to withstand the harsh conditions of sterilization. Our unique compilation of textiles is designed to be strong and soft for easy wrapping and durability. Manufactured for single use only. For maximum protection.

- Ideal for sterile storage and subsequently covering the hygienic work surface
- Unique composition of textile fibers makes the wrap soft enough for packaging yet robust for sterilisation and storage.




Safe for the environment, enhances safety for patients and users. IMS Monitor Tape - lead free and not made with natural rubber latex

- For labeling wrapped instrument cassettes
- The monitor tape is safe for the environment and excludes the risk of lead contamination
- Enhances safety for patients and users who suffer latex allergies
- Monitoring tape can be labeled as required
- Serves as a reliable process indicator: Black strips appear on the tape once the cassette has been sterilised in the autoclave.
- The IMS Monitor Tape - lead free and not made with natural rubber latex is simple to dispose of after use




WRAPPING TECHNIQUES

STANDARD GIFT



HOSPITAL STYLE



IMS STERILISATION WRAP

After the packed cassette has been opened, the sterilisation paper may be used as a base, thus forming a predominantly sterile environment. Additionally, contaminated cotton wool rolls, gauze and similar items can be disposed of hygienically in this paper.

Suitable for	Dimensions	Contents	Order No.
DIN Cassettes and Container Cassettes	600 x 600 mm	500 sheets	IMS-1210M
DIN Cassettes and Container Cassettes	600 x 600 mm	100 sheets	IMS-1210MH
1/4 DIN Cassettes	375 x 375 mm	1000 sheets	IMS-1211M
1/4 DIN Cassettes	375 x 375 mm	100 sheets	IMS-1211MH
1/2 DIN Cassettes	500 x 500 mm	500 sheets	IMS-1212M
1/2 DIN Cassettes	500 x 500 mm	100 sheets	IMS-1212MH

IMS BONDED DOUBLE WRAP

Suitable for	Dimensions	Contents	Order No.
1/4 DIN Cassettes	380 x 380 mm	480 pcs.	IMS-2215
1/2 DIN Cassettes	510 x 510 mm	480 pcs.	IMS-2220
DIN Cassettes and Container Cassettes	610 x 610 mm	240 pcs.	IMS-2224

IMS MONITOR TAPE

The IMS monitor tapes guarantee that the packaging does not open during the sterilisation process. Black strips appear on the tape once the cassette has been sterilised in the autoclave. The IMS Monitor Tape IMS-1460 - lead free and not made with natural rubber latex enhances safety for patients and users.

Lead and Natural Rubber Free Monitor Tape 3/4" (19 x 55 mm) [IMS-1460](#)

Also available: Monitor tape color-coded, blank or pre-printed (with specific procedural set-ups).

IMS MONITOR TAPE DISPENSER

Holds up to 9 rolls. Dispenser can also be mounted under cabinet.

up to 9 rolls, 27,3 x 15,9 x 19,4 cm [IMS-1270](#)
up to 4 rolls, 11,1 x 15,9 x 19,4 cm [IMS-1270H](#)

