

SINGLE-USE OPERATORY SUPPLIES




Used everyday in your operatory, single-use disposable products help prevent cross contamination and increase efficiency.

PATIENT TOWELS, BIBS, AND LINERS CROSSTEX

For over 60 years, we have been the benchmark of quality in patient towels, bibs, and liners. Our unique manufacturing process:

- Provides the often imitated, but never duplicated, Crosstex embossing pattern, which allows moisture to flow horizontally across the towel, protecting your patient's clothing.
- Offers our exclusive reinforced edge embossing design for added strength and durability where towel clips are placed.
- Bonds the layers of tissue (and poly, when applicable) securely together, eliminating separation of the layers.
- Creates the strength of the towel horizontally through the combination of our proprietary manufacturing equipment and tissue formation. This prevents tearing when a towel clip is applied.



SANI-TAB™ CHAIN-FREE PATIENT TOWEL

Don't let infection spread at a dangerous clip!

Dental bib clips and chains have been shown to pick up dangerous microorganisms from skin, saliva and plaque, presenting a potential source of cross-contamination in the dental operatory. Reduce the risk with Sani-Tab™ Chain-Free Patient Towels — the only chain-free patient towel with patented adhesive tabs eliminating the need for bib chains and clips. Sani-Tab maximizes infection control and patient comfort while saving time and money.

- Easy to use — simply position towel anywhere on the patient and apply pressure to the tab.
- Bib stays where placed to maximize coverage, unlike bib clips and other disposable towel holders.
- Can be placed high on patient to protect clothing.



TOWEL DISPENSER

- Made from polished acrylic.
- Sits on counter/shelf or easily mounts to wall with supplied mounting kit.
- Accommodates folded towels up to 19 in x 13 in (48.3 cm x 33.0 cm).
- Holds approximately 50 towels.
- Front opening allows for easier grasping of towels.



TOWEL DISPENSER

CODE	DETAILS	QTY
PTC	5 in x 6 in x 7 in (12.7 cm x 15.2 cm x 17.7 cm)	Each

PATIENT TOWELS, BIBS, AND LINERS CROSSTEX

3-Ply Tissue with 1-Ply Poly

DESCRIPTION	QTY	SIZE	FOLD	EMBOSSING PATTERN	CODE, COLOR		
Polyback™ Towels	500/Box	19 in x 13 in (48.3 cm x 33.0 cm)	Crosstex	Horizontal	WPXAQ , Aqua WPXBG , Beige WPXBL , Blue WPXDR , Dusty Rose	WPXGR , Green WPXLV , Lavender WPXPE , Peach WPXSI , Silver Grey	WPXWH , White WPXYE , Yellow WPXFP , Happy Feet WPXTB , Teddy Bear
Poly-Gard™ Towels	500/Box	19 in x 16 in (48.3 cm x 40.6 cm)	Crosstex	Horizontal	WGXBL , Blue WGXWH , White		
Poly-Gard™ Deluxe Towels	300/Box	18 in x 30 in (45.7 cm x 76.2 cm)	Vantage	Pebble	W3D1930BL , Blue W3D1930WH , White		
Polyback™ Contour Rib Towels	400/Box	18 in x 22 in (45.7 cm x 55.9 cm)	Vantage	Pebble	WERUECBL , Blue WERUECDR , Dusty Rose	WERUECGR , Green WERUEC , White	
Sani-Tab™ Chain-Free Patient Towel	400/Box	19 in x 13 in (48.3 cm x 33.0 cm)	Crosstex	Horizontal	WPXBLT , Blue WPXDRT , Dusty Rose WPXWHT , White		

2-Ply Tissue with 1-Ply Poly

DESCRIPTION	QTY	SIZE	FOLD	EMBOSSING PATTERN	CODE, COLOR		
Econoback™ Towels	500/Box	19 in x 13 in (48.3 cm x 33.0 cm)	Crosstex	Horizontal	WEXAQ , Aqua WEXBG , Beige WEXBL , Blue WEXDR , Dusty Rose WEXGR , Green	WEXLV , Lavender WEXPE , Peach WEXSI , Silver Grey WEXWH , White WEXYE , Yellow	WEXFP , Happy Feet WEXTB , Teddy Bear
Ultra-Gard™ Towels	500/Box	19 in x 16 in (48.3 cm x 40.6 cm)	Crosstex	Horizontal	WUXBL , Blue WUXYE , Yellow	WUXWH , White	
Sani-Tab™ Chain-Free Patient Towel	400/Box	19 in x 13 in (48.3 cm x 33.0 cm)	Crosstex	Horizontal	WEXBLT , Blue WEXWHT , White	WEXDRT , Dusty Rose	WEXLVT , Lavender

Extra Heavy Tissue with 1-Ply Poly

DESCRIPTION	QTY	SIZE	FOLD	EMBOSSING PATTERN	CODE, COLOR		
Proback™ Disposable Paper Towels	500/Box	19 in x 13 in (48.3 cm x 33.0 cm)	Crosstex	Horizontal	WEPBAQ , Aqua WEPBBL , Blue WEPBDR , Dusty Rose WEPBGR , Green	WEPBLV , Lavender WEPBPE , Peach WEPBSI , Silver Grey WEPBWH , White	WEPBYE , Yellow WEPBTB , Teddy Bear

3-Ply Tissue

DESCRIPTION	QTY	SIZE	FOLD	EMBOSSING PATTERN	CODE, COLOR		
Professional™ Regular Towels	500/Box	19 in x 13 in (48.3 cm x 33.0 cm)	Crosstex	Horizontal	WTXAQ , Aqua WTXBG , Beige WTXBL , Blue	WTXDR , Dusty Rose WTXGR , Green WTXWH , White	WTXYE , Yellow WTXTB , Happy Feet WTXSI , Silver Grey

SALIVA EJECTORS CROSSTEX

COMFORT PLUS™ PREMIUM SALIVA EJECTORS

- The kinder, gentler saliva ejector...takes in fluid, not tissue!
- Patented tip design enhances patient comfort – will not grab mucosal tissue.
- Facilitates uninterrupted fluid and light debris removal.
- Optimizes fluid removal as patented tip design prevents tissue from blocking suction.
- Holds shape when formed as result of proprietary manufacturing process.
- Available in multiple color combinations, with and without scent.



COMFORT PLUS™ PREMIUM SALIVA EJECTORS

CODE	DETAILS	QTY
ZWWCP	White/White	100/Bag
ZCWCP	Clear/White	100/Bag
ZWPCPBG	White/Pink Bubblegum scented	100/Bag
ZWGCPM	White/Green Mint scented	100/Bag

HYGOFORMIC™ SALIVA EJECTORS

- Combination saliva ejector and tongue retractor.
- Unique coil design is adjustable for maximum patient comfort.
- Smooth edges eliminate irritation – ideal for patients with sensitive tissue.
- Environmentally friendly – biodegradable.



Made outside of the USA.

HYGOFORMIC™ SALIVA EJECTORS

CODE	DETAILS	QTY
ZFF	White	100/Bag
ZFFAD	Adapters	10/Bag

PREMIUM™ SALIVA EJECTORS AND ADVANTAGE™ SALIVA EJECTORS

- Holds shape when formed as result of proprietary manufacturing process.
- Soft, non-removable tips are bonded to the tube for maximum patient safety.
- Provides optimal suction.
- Matte finish eliminates glare.
- Available in multiple color combinations.
- Advantage Saliva Ejectors have identical manufacturing process to Premium Saliva Ejectors except with a glossy finish.



PREMIUM™ SALIVA EJECTORS

CODE	DETAILS	QTY
ZCB	Frosted/Blue	100/Bag
ZCW	Frosted/Blue	100/Bag
ZWW	White/White	100/Bag
ZPP	Pink/Pink	100/Bag
ZCMC	Multi Color	100/Bag

ADVANTAGE™ SALIVA EJECTORS

CODE	DETAILS	QTY
ZBBI	Blue/Blue	100/Bag
ZCBI	Clear/Blue	100/Bag
ZCCI	Clear/Clear	100/Bag
ZWWI	White/White	100/Bag

EVACUATOR TIPS CROSSTEX

HVE COMBO-TIP™ EVACUATORS

- Manufactured from high density polyethylene.
- 2 tips in one – vented on one end/non-vented on the other.
- Unique “S” tip design eliminates sharp edges enhancing patient safety and comfort.
- Stiff, firm construction allows for cheek retraction.



HVE COMBO-TIP™ EVACUATORS

CODE	DETAILS	QTY
ZET	Combo Tip™, White	100/Bag

MAXVAC PLUS WITH COLOR™ COMBO TIP™

- Identical to the Combo-Tip™ in design but not in color.
- Available in 5 décor-pleasing colors.



MAXVAC PLUS™ COMBO TIP™ EVACUATORS

CODE	DETAILS	QTY
ZETMVBL	Blue	100/Bag
ZETMVGR	Green	100/Bag
ZETMVPK	Pink	100/Bag
ZETMVLV	Purple	100/Bag
ZETMVYE	Yellow	100/Bag
ZETMVP	Assorted (2 bags/color per case)	10 Bags/ Case

NON-VENTED EVACUATOR TIPS

- Identical to Combo-Tip™ Evacuators except both ends are non-vented.



NON-VENTED EVACUATOR TIPS


CODE	DETAILS	QTY
ZETNV	Non-vented, White	50/Bag

EVACUATOR TIPS CROSSTEX

SPARKLE™ AIR WATER SYRINGE TIPS

- Requires no adapter – locking groove ensures tip is securely installed.
- Quad Core™ Air Channel ensures complete air/water separation.
- Translucent colors allow for visualization of water flow/blockage.
- Beveled edge provides increased patient comfort and safety.
- Rounded non-scratch end prevents O-ring damage.
- Single-use eliminates risk of cross-contamination and cleaning of syringe between patients.
- Quick, easy, economical.
- Available in 5 translucent colors and opaque white.

Made outside of the USA.




SPARKLE™ AIR WATER SYRINGE TIPS

CODE	DETAILS	QTY
BCSAWSCL	Clear	250/Bag
BCSAWSWH	Green	250/Bag
BCSAWS	Assorted (Blue, Purple, Yellow, Amber)	250/Bag
BCSAWSB	Assorted (Blue, Purple, Yellow, Amber)	1500/Bag

SURGICAL ASPIRATOR TIPS

- Manufactured from durable, high quality polyethylene with a 30° angle.
- Tips are designed to fit all high-volume evacuation lines.
- Strong construction enables tip to be used for cheek retraction as it will not bend.
- Available in 3 tip widths.

Made outside of the USA.



SURGICAL ASPIRATOR TIPS

CODE	DETAILS	QTY
ZSABL	Small Tip, Blue 0.0625 in (0.16 cm)	25/Bag
ZSAWH	Standard Tip, White 0.125 in (0.32 cm)	25/Bag
ZSAGR	Large Tip, Green 0.25 in (0.64 cm)	25/Bag

COTTON AND NON-WOVENS CROSSTEX

COTTON ROLLS

- Made from 100% cotton.
- Contains no cellulose or rayon fibers.
- Soft, pliable, non-linting.
- Will not adhere to mucous membrane.
- Maintains shape when wet.
- Absorbs up to 9 times their weight.
- Available in Premium and Advantage Plus.
- Non-sterile.

Made outside of the USA.



PREMIUM™ COTTON ROLLS

CODE	DETAILS	QTY
DNC	#2 Medium 1.5 in x 0.375 in (3.8 cm x 1 cm)	2,000/ Box

ADVANTAGE PLUS™ COTTON ROLLS

CODE	DETAILS	QTY
DNAC	#2 Medium 1.5 in x 0.375 in (3.8 cm x 1 cm)	2,000/ Box

COTTON ROLLS, BRAIDED NON-STERILE

- Made from 100%, short fiber cotton, non-sterile.
- Bound by thin strands of polyester filament.
- Starch-free for maximum softness.
- Tightly braided to maintain shape when wet.
- Highly absorbent.
- Will not adhere to mucous membrane.

Made outside of the USA.



COTTON ROLLS, BRAIDED NON-STERILE

CODE	DETAILS	QTY
DNCB	1.5 in x 0.375 in (3.8 cm x 1 cm)	2,000/ Box

COTTON ROLL DISPENSER DELUXE

- Flip-up top for easy loading.
- Holds up to 50 cotton rolls.
- Slide out drawer dispenses up to 4 rolls.
- Non-skid base for easy dispensing.



COTTON ROLLS, BRAIDED NON-STERILE

CODE	DETAILS	QTY
PDCBL	Blue 4 in x 3.3 in x 2 in (10.2 cm x 8.4 cm x 5.1 cm)	1 Each
PDCWH	White 4 in x 3.3 in x 2 in (10.2 cm x 8.4 cm x 5.1 cm)	1 Each

COTTON AND NON-WOVENS CROSSTEX

COTTON TIPPED APPLICATORS

- 100% cotton tip.
- All wood dowels.
- Extremely tight bud minimizes linting.
- Applicators are sealed in autoclavable bags of 100.
- Available in Premium™ Cotton Tipped Applicators (domestic) and Advantage Plus™ Cotton Tipped Applicators (imported).



PREMIUM™ COTTON TIPPED APPLICATORS

CODE	DETAILS	QTY
H3C	3 in (7.6 cm)	1,000/Box
H6C	6 in (15.2 cm)	1,000/Box

ADVANTAGE PLUS™ COTTON TIPPED APPLICATORS

CODE	DETAILS	QTY
H3AC	3 in (7.6 cm)	1,000/Box
H6AC	6 in (15.2 cm)	1,000/Box

Made outside of the USA.

PREMIUM™ ALL-GAUZE SPONGES

- Made from USP 100% gauze.
- Available in 8-ply, 12-ply or 16-ply.
- Highly absorbent.
- Non-sterile.

Made outside of the USA.



PREMIUM™ NON-STERILE ALL-GAUZE SPONGES

CODE	DETAILS	QTY
ENC212	12-Ply 2 in x 2 in (5.1 cm x 5.1 cm)	8,000/Case (200/Sleeve)
ENC2	8-Ply 2 in x 2 in (5.1 cm x 5.1 cm)	5,000/Case (200/Sleeve)
ENC312	12-Ply 3 in x 3 in (7.6 cm x 7.6 cm)	4,000/Case (200/Sleeve)
ENC38	8-Ply 3 in x 3 in (7.6 cm x 7.6 cm)	4,000/Case (200/Sleeve)
ENC416	16-Ply 4 in x 4 in (10.2 cm x 10.2 cm)	2,000/Case (200/Sleeve)
ENC412	12-Ply 4 in x 4 in (10.2 x 10.2 cm)	2,000/Case (200/Sleeve)
ENC48	8-Ply 4 in x 4 in (10.2 x 10.2 cm)	4,000/Case (200/Sleeve)

ULTRA-GAUZE™ NON-WOVEN SPONGES

- 4-ply, non-sterile.
- 40 g weight (weight per m2 of rayon/poly fabric).
- More absorbent than traditional sponges, exclusive 70/30 polyester formed non-woven material replicates the look of woven gauze.
- Ideal for placement on incision/extraction sites as it easily releases, for increased patient comfort.
- Extremely soft, smooth, strong—virtually lint-free.
- Ideal for sustained patient use, general cleaning.



ULTRA-GAUZE™ NON-WOVEN SPONGES

CODE	DETAILS	QTY
ENCNWU	2 in x 2 in (5.1 cm x 5.1 cm)	5,000/Case (200/Sleeve)
ENCNWUB	2 in x 2 in (5.1 cm x 5.1 cm)	8,000/Case (200/Sleeve)
ENC3NWU	3 in x 3 in (7.6 cm x 7.6 cm)	4,000/Case (200/Sleeve)
ENC4NWU	4 in x 4 in (10.2 cm x 10.2 cm)	2,000/Case (200/Sleeve)

COTTON AND NON-WOVENS CROSSTEX

COTTON-FILLED SPONGES

- Made from USP cotton gauze.
- Pure, virgin white cotton with 100% cotton-fiber filling.
- Highly absorbent.
- Available in sterile and non-sterile.

Made outside of the USA.



PREMIUM™ COTTON-FILLED SPONGES, EXODONTIA, NON-STERILE

CODE	DETAILS	QTY
ENCC	2 in x 2 in (5.1 cm x 5.1 cm) 125 g weight	20,000/Case (200/Sleeve)
ENCC3	3 in x 3 in (7.6 cm x 7.6 cm) 125 g weight	4,000/Case (200/Sleeve)
ENCC4	4 in x 4 in (10.2 cm x 10.2 cm) 125 g weight	2,000/Case (200/Sleeve)

PREMIUM™ COTTON-FILLED SPONGES, EXODONTIA, STERILE

CODE	DETAILS	QTY
ESCC	2 in x 2 in (5.1 cm x 5.1 cm) 125 g weight	20,000/Case (200/Sleeve)

PREMIUM™ COTTON-FILLED SPONGES, EXODONTIA, STERILE

CODE	DETAILS	QTY
ENACC	2 in x 2 in (5.1 cm x 5.1 cm) 90 g weight	20,000/Case (200/Sleeve)

NON-WOVEN SPONGES

- 4 Ply, Non-Sterile.
- Unique rayon/polyester blend.
- Provides maximum absorbency.
- Less adhesion to wounds for faster healing.
- Soft, smooth, strong, virtually lint-free.
- Ideal for sustained patient use, disinfecting, general cleaning.
- Available in Premium™ Non-Wovens, Advantage Plus™ Non-Wovens, and Advantage™ Non-Wovens in multiple gram weights and sizes.



PREMIUM™ NON-WOVENS, NON-STERILE

CODE	DETAILS	QTY
ENCNW4K	2 in x 2 in (5.1 cm x 5.1 cm) 50 g weight	4,000/Case (200/Sleeve)
ENCNW	2 in x 2 in (5.1 cm x 5.1 cm) 50 g weight	20,000/Case (200/Sleeve)
ENC3NW	3 in x 3 in (7.6 cm x 7.6 cm) 40 g weight	20,000/Case (200/Sleeve)
ENC4NW	4 in x 4 in (10.2 cm x 10.2 cm) 40 g weight	20,000/Case (200/Sleeve)

ADVANTAGE PLUS™ NON-WOVENS, NON-STERILE

CODE	DETAILS	QTY
ENCNWL	2 in x 2 in (5.1 cm x 5.1 cm) 35 g weight	20,000/Case (200/Sleeve)

ADVANTAGE PLUS™ NON-WOVENS, NON-STERILE

CODE	DETAILS	QTY
ENCNWL A	2 in x 2 in (5.1 cm x 5.1 cm) 30 g weight	5,000/Case (200/Sleeve)
ENC3NWL A	3 in x 3 in (7.6 cm x 7.6 cm) 30 g weight	4,000/Case (200/Sleeve)
ENC4NWL A	4 in x 4 in (10.2 cm x 10.2 cm) 30 gram weight	2,000/Case (200/Sleeve)

COTTON AND NON-WOVENS CROSSTEX

ULTRA-PURE™ NON-WOVEN SPONGES

- 4-ply, non-sterile.
- Biodegradable – a great way to reduce your carbon footprint!
- Made from 100% USA-grown cotton, providing increased comfort and faster absorbency when compared to rayon/ polyester and cotton blends.
- Hypoallergenic, no toxic chemicals used in manufacturing process; all natural whitening process (H2O2 whitening) without use of bleach and/or chlorine.
- Exclusive automated manufacturing process reduces cross-contamination risks.
- Ideal for placement on incision/extraction sites as it easily releases, for increased patient comfort.
- Lint-free – extremely soft, smooth, strong and will not adhere to instruments.

Made outside of the USA.



ULTRA-PURE™ NON-WOVEN SPONGES

CODE	DETAILS	QTY
ENC100	2 in x 2 in (5.1 cm x 5.1 cm) 45 g weight	200/Sleeve
ENC4100	4 in x 4 in (10.2 cm x 10.2 cm) 40 g weight	200/Sleeve

Cotton & Non-Woven Sponge Chart

PRODUCT	CODE	STYLE	PLY	STERILITY	SIZE	WEIGHT	QTY
Premium™ Gauze Sponges	ENC212	100% Gauze	12-Ply	Non-Sterile	2 in x 2 in (5.1 cm x 5.1 cm)	-	8,000/Case
	EN212	100% Gauze	8-Ply	Non-Sterile	2 in x 2 in (5.1 cm x 5.1 cm)	-	5,000/Case
	ENC312	100% Gauze	12-Ply	Non-Sterile	3 in x 3 in (7.6 cm x 7.6 cm)	-	4,000/Case
	ENC38	100% Gauze	8-Ply	Non-Sterile	3 in x 3 in (7.6 cm x 7.6 cm)	-	4,000/Case
	ENC416	100% Gauze	16-Ply	Non-Sterile	4 in x 4 in (10.2 cm x 10.2 cm)	-	2,000/Case
	ENC412	100% Gauze	12-Ply	Non-Sterile	4 in x 4 in (10.2 cm x 10.2 cm)	-	2,000/Case
	ENC48	100% Gauze	8-Ply	Non-Sterile	4 in x 4 in (10.2 cm x 10.2 cm)	-	4,000/Case
Premium™ Non-Sterile Sponges	ENCC	Cotton-filled	-	Non-Sterile	2 in x 2 in (5.1 cm x 5.1 cm)	125 g	5,000/Case
	ENCC3	Cotton-filled	-	Non-Sterile	3 in x 3 in (7.6 cm x 7.6 cm)	125 g	4,000/Case
	ENCC4	Cotton-filled	-	Non-Sterile	4 in x 4 in (10.2 cm x 10.2 cm)	125 g	2,000/Case
Premium™ Sterile Sponges	ESCC	Cotton-filled	-	Sterile	2 in x 2 in (5.1 cm x 5.1 cm)	125 g	5,000/Case
Advantage Plus™ Sponges	ENACC	Cotton-filled	-	Non-Sterile	2 in x 2 in (5.1 cm x 5.1 cm)	90 g	5,000/Case
Premium™ Sponges	ENCNW	Non-Woven	4-Ply	Non-Sterile	2 in x 2 in (5.1 cm x 5.1 cm)	50 g	5,000/Case
	ENC3NW	Non-Woven	4-Ply	Non-Sterile	3 in x 3 in (7.6 cm x 7.6 cm)	40 g	4,000/Case
	ENCNW4K	Non-Woven	4-Ply	Non-Sterile	2 in x 2 in (5.1 cm x 5.1 cm)	50 g	4,000/Case
	ENC4NW	Non-Woven	4-Ply	Non-Sterile	4 in x 4 in (10.2 cm x 10.2 cm)	40 g	2,000/Case
Advantage Plus™ Sponges	ENCNWL	Non-Woven	4-Ply	Non-Sterile	2 in x 2 in (5.1 cm x 5.1 cm)	35 g	5,000/Case
Advantage™ Sponges	ENCNWLA	Non-Woven	4-Ply	Non-Sterile	2 in x 2 in (5.1 cm x 5.1 cm)	30 g	5,000/Case
	ENC3NWLA	Non-Woven	4-Ply	Non-Sterile	3 in x 3 in (7.6 cm x 7.6 cm)	30 g	4,000/Case
	ENC4NWLA	Non-Woven	4-Ply	Non-Sterile	4 in x 4 in (10.2 cm x 10.2 cm)	30 g	2,000/Case
Ultra-Gauze™ Sponges	ENCNWU	Non-Woven	4-Ply	Non-Sterile	2 in x 2 in (5.1 cm x 5.1 cm)	40 g	5,000/Case
	ENCNWUB	Non-Woven	4-Ply	Non-Sterile	2 in x 2 in (5.1 cm x 5.1 cm)	40 g	8,000/Case
	ENC3NWU	Non-Woven	4-Ply	Non-Sterile	3 in x 3 in (7.6 cm x 7.6 cm)	40 g	4,000/Case
	ENC4NWU	Non-Woven	4-Ply	Non-Sterile	4 in x 4 in (10.2 cm x 10.2 cm)	40 g	2,000/Case
Ultra-Pure™ Sponges	ENC100	Non-Woven	4-Ply	Non-Sterile	2 in x 2 in (5.1 cm x 5.1 cm)	45 g	5,000/Case
	ENC4100	Non-Woven	4-Ply	Non-Sterile	4 in x 4 in (10.2 cm x 10.2 cm)	40 g	2,000/Case

OPERATORY PROTECTION_{CROSSTEX}

BRACKET TRAY COVERS


- Manufactured from a unique, heavy weight paper.
- Impedes flow of moisture, minimizing seep-through.
- Available in multiple sizes and colors.

BRACKET TRAY COVERS, LIQUID-PROOF

- Polyethylene-backed paper prevents seep-through.
- Available in two sizes, white only.

DUAL USE PRODUCT:

- Poly-side down – provides positive moisture barrier.
- Poly-side up – provides ideal surface for quick mixing.



Bracket Tray Covers

DESCRIPTION	QTY	SIZE	CODE, COLOR
Weber/Chayes Bracket Tray Covers	1000/Carton	A 9.5 in x 12.25 in (24.1 cm x 31.1 cm)	FABL , Blue FAWH , White FADR , Dusty Rose FAGR , Green
Ritter Bracket Tray Covers	1000/Carton	B 8.5 in x 12.25 in (21.6 cm x 31.1 cm)	FBBL , Blue FBWH , White FBDR , Dusty Rose FBGR , Green FBAQ , Aqua FBPE , Peach FBYE , Yellow FBBG , Beige FBSI , Silver Grey FBLV , Lavender
Weber Hi-Lo Bracket Tray Covers	1000/Carton	C 11 in x 17.5 in (27.9 cm x 44.5 cm)	FCBL , Blue FCWH , White FCGR , Green
SS White #20 Bracket Tray Covers	1000/Carton	D 10.25 in x 15.75 in (26 cm x 40 cm)	FDBL , Blue FDWH , White FDGR , Green
Midwest Bracket Tray Covers	1000/Carton	E 9 in x 13.5 in (22.9 cm x 34.3 cm)	FEBL , Blue FEWH , White
Mini Bracket Tray Covers	1000/Carton	M 5 in x 8 in (12.7 cm x 20.2 cm)	FMINI , White
Standard Round Bracket Tray Covers	1000/Carton	O 13.5 in diameter (34.3 cm diameter)	FOBL , Blue FOWH , White
Viking Bracket Tray Covers	1000/Carton	V 8.5 in x 11 in (21.6 cm x 27.9 cm)	FVWH , White
Cox Bracket Tray Covers	1000/Carton	X 11 in x 11 in (27.9 cm x 27.9 cm)	FCOX , White

Bracket Tray Covers, Liquid Proof


DESCRIPTION	QTY	SIZE	CODE, COLOR
Ritter Liquid Proof Bracket Tray Covers	500/Carton	B 8.5 in x 12.25 in (21.6 cm x 31 cm)	FBC , White
Adec Liquid Proof Bracket Tray Covers	500/Carton	F 8.5 in x 11 in (21.6 cm x 27.9 cm)	FF , White

OPERATORY PROTECTIONCROSSTEX

HEADREST COVERS, NON-WOVEN

- Soft, breathable non-woven fabric provides reliable protection against cross-contamination.
- Provides maximum patient comfort as fabric is quiet – does not rustle.
- Available in 2 sizes.

Made outside of the USA.



STANDARD HEADREST COVERS, NON-WOVEN 10 IN X 10 IN (25.4 CM X 25.4 CM)

CODE	DETAILS	QTY
LOCB	Blue	500/Case
LOC	White	500/Case
LOCL	Lavender	500/Case

EXTRA LARGE HEADREST COVERS, NON-WOVEN 10 IN X 13 IN (25.4 CM X 33 CM)

CODE	DETAILS	QTY
L3CB	Blue	500/Case
L3C	White	500/Case

HEADREST COVERS, PLASTIC

- Clear, soft 1 mil thick polyethylene provides impervious barrier between the patient and the headrest for maximum protection from cross-contamination.
- 3/8 inch (0.95 cm) gusset design allows cover to lay flat on headrest for added patient comfort.
- Available in two sizes.



REGULAR HEADREST COVERS, PLASTIC 9.5 IN X 11 IN (24.1 CM X 27.9 CM)

CODE	DETAILS	QTY
LOCP	Clear	250/Dispenser Box
LOCPW	White	250/Dispenser Box

JUMBO HEADREST COVERS, PLASTIC 9.5 IN X 14 IN (24.1 CM X 35.6 CM)

CODE	DETAILS	QTY
L3CP	Clear	250/Dispenser Box
L3CPW	White	250/Dispenser Box

HEADREST COVERS, POLYCOATED

- Soft facial tissue combined with poly to provide positive fluid barrier for maximum protection from cross-contamination.
- Soft, strong and absorbent.
- Available in multiple sizes and colors.



DESCRIPTION	QTY	SIZE	CODE, COLOR		
Standard Headrest Covers, Polycoated	500/Case	10 in x 10 in (25.4 cm x 25.4 cm)	LOBL, Blue LOWH, White	LODR, Dusty Rose LOGR, Green	LOLV, Lavender
Extra Large Headrest Covers, Polycoated	500/Case	10 in x 13.0 in (25.4 cm x 33 cm)	L3BL, Blue L3WH, White	L3DR, Dusty Rose L3GR, Green	L3LV, Lavender
Jumbo Headrest Covers, Polycoated	500/Case	13 in x 13 in (33.0 cm x 33.0 cm)	L13BL, Blue L13WH, White		
Extra Jumbo Headrest Covers, Polycoated	500/Case	11.5 in x 15 in (29.2 cm x 38.1 cm)	L15WH, White		

OPERATORY PROTECTIONCROSSTEX

BARRIER FILM WITH FINGER-LIFT EDGE

- Unique 3/16 in (0.476 cm) non-adhesive finger-lift edge allows quick, easy placement/removal.
- Light tack adhesive easily attaches to surfaces without leaving residue.
- Ideal for areas that are hard-to-reach, clean or disinfect.
- Available in blue and clear.



BARRIER FILM WITH FINGER-LIFT EDGE

CODE	DETAILS	QTY
BFBL	4 in x 6 in (10.2 cm x 15.2 cm), Blue	1200 Sheets/Roll
BFCL	4 in x 6 in (10.2 cm x 15.2 cm), Clear	1200 Sheets/Roll

PILLOW CASES

- Soft facial tissue combined with poly to provide positive fluid barrier for maximum protection from cross-contamination.
- Soft, strong and absorbent.

Made outside of the USA.



PILLOW CASES

CODE	DETAILS	QTY
PC2130	21 in x 30 in (53.3 cm x 76.2 cm)	100 Sleeves/Box

CHAIR SLEEVES

- Universal, clear plastic sleeves for patient chairs.
- Half sleeve design covers chair back and side chair controls.
- Reduces risk of cross-contamination for airborne aerosols and body fluids.
- Quick, easy, economical.

Made outside of the USA.



CHAIR SLEEVES

CODE	DETAILS	QTY
BCCS	Clear 27.5 in x 24 in (69.9 cm x 61 cm)	225 Sleeves/Roll

OPERATORY PROTECTIONCROSSTEX

LIGHT SLEEVES, T-HANDLE

- Patented clear plastic sleeve for T-style light handles.
- White indicator arrow for quick identification of die-cut location.
- Will not come off during use.
- Quick, easy, economical.

Made outside of the USA.



LIGHT SLEEVES, T-HANDLE

CODE	DETAILS	QTY
BCLST	Clear with white arrow 4 in x 5.75 in (10.2 x 14.6 cm)	500/Box

X-RAY SLEEVE

- Universal clear plastic sleeve fits most x-ray heads.
- Reduces risk of cross-contamination for airborne aerosols and body fluids.
- Will not adversely affect quality of the x-ray image.
- Quick, easy, economical.

Made outside of the USA.



X-RAY SLEEVE

CODE	DETAILS	QTY
BCXR	Clear 15 in x 26 in (38.1 cm x 66 cm)	250/Box

SYRINGE SLEEVES

- Clear plastic sleeve easily fits over 3-way syringe, high/low-speed hand pieces and/or air abrasion hand pieces.
- Reduces risks of cross-contamination for airborne aerosols and body fluids.

Available in three configurations:

- Without opening – simply place over air/ water syringe or handpiece and push through.
- With opening – convenient pre-cut opening for easy insertion.
- With opening and adhesive – Convenient, pre-cut opening for easy insertion with clear adhesive strip that attaches to tubing to prevent sliding of sleeve, creating a closed barrier.



SYRINGE SLEEVE WITHOUT OPENING

CODE	DETAILS	QTY
BCSS	Clear 2.5 in x 10 in (6.4 cm x 25.4 cm)	500/Box

SYRINGE SLEEVE WITH OPENING

CODE	DETAILS	QTY
BCSSO	Clear 2.5 in x 10 in (6.4 cm x 25.4 cm)	500/Box
BCSSBO	Blue 2.5 in x 10 in (6.4 cm x 25.4 cm)	500/Box

SYRINGE SLEEVE WITH OPENING AND CLEAR ADHESIVE STRIP

CODE	DETAILS	QTY
BCSSA	Clear 2.5 in x 11.4 in (6.4 cm x 29 cm)	500/Box

Made outside of the USA.

OPERATORY PROTECTION CROSSTEX

TRAY SLEEVES

- Clear plastic sleeve for patient trays.
- Secure fit end-fold provides maximum protection and reduces the risk of cross-contamination from airborne aerosols and body fluids.
- Quick, easy, economical.
- Available in 4 convenient sizes.



TRAY SLEEVES, ADEC/SIZE B-RITTER

CODE	DETAILS	QTY
BCTS1014	Clear 10.5 in x 14 in (26.7 cm x 35.6 cm)	500/Box

TRAY SLEEVES, MINI TRAYS

CODE	DETAILS	QTY
BCTS710	Clear 7.5 in x 10.5 in (19.1 cm x 26.7 cm)	500/Box

Made outside of the USA.

TRAY SLEEVES, CHAYES/SIZE B RITTER

CODE	DETAILS	QTY
BCTS1114	Clear 11.625 in x 14.5 in (29.5 cm x 36.8 cm)	500/Box

Made outside of the USA.

TRAY SLEEVES, MIDWEST/WEBER

CODE	DETAILS	QTY
BCTS1116	Clear 11.5 in x 16 in (29.2 cm x 40.6 cm)	500/Box

Made outside of the USA.

LAB TRANSPORT BAGS

- Meets OSHA requirements and recommendation of NCCLS (National Committee for Clinical Laboratory Standards) for transport of medical specimens and/or dental lab cases.
- Convenient two-pouch system keeps Rx/lab slip separate from specimens/lab case as required by OSHA.
- Biohazard symbol for easy identification and OSHA compliance.
- Permanent adhesive closure does not contact gloved hand – tamper-evident.
- White front can be written on, allowing for easy recording of information.
- Quick, easy, economical.

Made outside of the USA.



LAB TRANSPORT BAGS

CODE	DETAILS	QTY
BAGS	6 in x 10 in (15.2 cm x 25.4 cm)	250/Box

CUPS CROSSTEX

PLASTIC CUPS

- Unique embossing from greater strength and easy dispensing.
- Double-rolled lip for easy separation and enhanced patient comfort.
- Thick-walled for added strength and patient safety.
- Double-coated for lustrous, high-gloss finish.
- Available in 9 colors.



PLASTIC CUPS 5 OZ (148 ML)

CODE	DETAILS	QTY
CXBG	Beige	1000/Case
CXBL	Blue	1000/Case
CXDR	Dusty Rose	1000/Case
CXGR	Green	1000/Case
CXLV	Lavender	1000/Case
CXSI	Silver Grey	1000/Case
CXWH	White	1000/Case
CXYE	Yellow	1000/Case
CXCL	Clear/Translucent	2500/Case

CUP DISPENSER

- Contemporary design fits décor of every office.
- Dispenses fifty 5-ounce (148 ml) plastic and wax-coated paper cups—one cup at a time.
- Includes self-adhesive strips for easy wall mounting.

Made outside of the USA.



CUP DISPENSER

CODE	DETAILS	QTY
PCC	White base with clear top	1 Each

MEDICINE / MIXING CUPS

- Features fluid ounce, dram, cc, ml, teaspoon and tablespoon units with clearly marked graduations.
- Convenient for dispensing or mixing.
- Unbreakable polypropylene material can be squeezed to form a pour spout.



MEDICINE/MIXING CUPS

CODE	DETAILS	QTY
CX1	1 fl oz (30 ml)	100/Sleeve

SURGICAL LINE GLOVES OMNIA

Anatomically-designed, sterile surgical gloves made of latex with a low content of allergenic proteins and presence of powder.

CODE	DETAILS	MATERIAL and SIZE	PAIRS
32.G5000.00	Pre powdered surgical gloves - Digitosteril	latex - size 6	50
32.G5001.00	Pre powdered surgical gloves - Digitosteril	latex - size 6 ^{1/2}	50
32.G5002.00	Pre powdered surgical gloves - Digitosteril	latex - size 7	50
32.G5003.00	Pre powdered surgical gloves - Digitosteril	latex - size 7 ^{1/2}	50
32.G5004.00	Pre powdered surgical gloves - Digitosteril	latex - size 8	50
32.G5005.00	Pre powdered surgical gloves - Digitosteril	latex - size 8 ^{1/2}	50
32.G5020.00	Sterile surgical gloves without powder	latex - size 6	50
32.G5021.00	Sterile surgical gloves without powder	latex - size 6 ^{1/2}	50
32.G5022.00	Sterile surgical gloves without powder	latex - size 7	50
32.G5023.00	Sterile surgical gloves without powder	latex - size 7 ^{1/2}	50
32.G5024.00	Sterile surgical gloves without powder	latex - size 8	50
32.G5025.00	Sterile surgical gloves without powder	latex - size 8 ^{1/2}	50



BIBS OMNIA

Salvet wipes are extra-absorbent, available in different colors: green, light blue and blue

CODE	DETAILS	COLOR
30.B0026.00	Salvet bib cm 32.5x45 cm	green
30.B0027.00	Salvet bib cm 32.5x45 cm	blue
30.B0022.00	Salvet bib cm 32.5x45 cm	dark blue
30.B0025.00	Salvet bib cm 32.2 x45 cm	orange
30.B0030.00	Salvet bib cm 32.2 x45 cm	lilac
30.B0023.00	Salvet bib cm 32.2 x45 cm	black
30.B0050.00	Salvet bib cm 32.2 x45 cm	teddy bear
30.B0028.00	Salvet bib cm 32.2 x45 cm	yellow



OMNISORB COTTON ROLLS OMNIA



The special function of these 100% pure cotton high-absorption rolls is to separate the cheeks from the teeth to prepare an easy access to the operating field. These rolls must not be used to absorb the saliva. Omnia also offers a practical cotton roll dispenser with a drawer.

CODE	DETAILS	QTY
30.Z1008.00	OmniSorb 8 mm Ø pure cotton rolls	300 g.
30.Z1009.00	OmniSorb 10 mm Ø pure cotton rolls	300 g.
30.Z1007.00	OmniSorb 12 mm Ø pure cotton rolls	300 g.
30.Z1031.00	Cotton roll with automatic dispenser	1 piece

For more information on Hu-Friedy products or to place an order, call at 00800 48 37 43 39 or visit HuFriedyGroup.eu

For more information on Crosstex products, call +39 05 24 52 74 53 or visit Crosstex.com

For more information on Omnia products or to place an order, call +39 05 24 52 74 53 or visit OmniaSpa.eu

Hu-Friedy Mfg. Co., LLC • European Headquarters • Lyoner Str. 9 • 60528 Frankfurt am Main, Germany • HuFriedyGroup.eu

As some of the products included in this catalogue might not be available in some countries, we invite you to verify the availability of the product with your local authorised dealer or to contact our customer care at the free phone number 00800 48 37 43 39. Available while stocks last. All other company and product names are trademarks of Hu-Friedy Mfg. Co. LLC, its affiliates or related companies, unless otherwise noted.

©2021 Hu-Friedy Mfg. Co., LLC All rights reserved. 90.L4005.EN/0921