


OMNIA PROCEDURAL SETS

SIMPLIFY YOUR
DENTAL ROUTINE



MAKE THE RIGHT CHOICE

Alongside the spread of more and more diseases and infections in the world, there is an increase in the awareness of all **infection control** processes that have to be followed during **surgical** and **non-surgical procedures**.

We have the solutions that help you **minimize the exposure to cross-infection** and blood borne pathogens during daily procedures: choose from a selection of **single-use disposables** (masks, gowns, drapes, accessories, medical adhesives, gauzes, pads) which are assembled in **Omnia Sets**. You will **optimize time, organization and efficiency** during chairside preparation while increasing staff and patients' **safety**.

Omnia Sets consist of individually configured single-use items. For example:

- Surgical drapes
- Surgical gowns
- Irrigation and aspiration systems
- Protection of the dental unit
- Accessories for surgery (sheaths and protection for instruments and physiodispensers, medical adhesives, gauzes, pads and irrigation lines)

TYPES OF OMNIA PROCEDURAL SETS

NON STERILE

Necessary protection during all kinds of **general procedures**.



STERILE

Reliable protection during **routine and advanced surgery**.



CUSTOMIZED

Perfect configuration for **your specific needs**.



OMNIA PROCEDURAL SETS: THE BENEFITS



HIGH QUALITY

- Production of packs occurs in a controlled environment classified as “**ISO CLASS 8**” pursuant to EN ISO 14644-1.
- The packs are packaged in double sterile packaging, medical paper and blister, in extremely rigorous procedures to meet the **highest control standards**.
- Constant after-sales surveillance and vigilance have been implemented in order to ensure that the packs are always at the forefront in terms of safety, suitability for use and customer satisfaction.
- During the internal production phases, monitoring is direct and constant so as to ensure ever-increasing quality standards and safety.
- **Sterility guaranteed for 5 years** from the date of production (if the package is not damaged).



COMPLIANCE

- Sterile Packs quality is well-documented by the knowledge and **compliance with current regulations** relating to medical devices
- Procedural Sets are compliant to the stricter national and local regulations.
- The label for sterile packs contains all the informations to identify the device (code, lot, expiration date). There are two extra labels that allow **easy data storage** and **traceability** of all materials and components used.



EFFICIENCY

- Practical, compact and **ready to use**.
- Contains all the disposable items which the dental team needs for the different types of operations.
- All products are already stacked in the order they will be used: once inside the operating theatre the dental team will just have to unwrap the pack, thereby making the set-up time easier and **more efficient**, thus raising the level of professionalism of the performance.
- The time saved will allow more focus towards the clinical case and will also increase the number of procedures possible.



SAVINGS

- Omnia Sets can significantly reduce the time of the various activities which enable carrying out procedures.
- Procedural packs make **planning processes, material purchasing, administration, set-up** and **traceability faster**.
- Significant reduction of total costs for each procedure.
- **Cost-effective**.



SUSTAINABILITY

- Procedural packs dramatically **reduce the amount of packaging material** compared with individually packaged products.
- Despite their use producing waste every day, they are preferable to reusable materials. The safety of single use items is guaranteed at every use and they don't have to be repeatedly washed and sterilised by medical personnel, with consequent saving of resources.

NON STERILE

OMNIA PROTECTION SETS

Protection Sets for Operator, Assistant and Patient during all kinds of general procedures.



Premium High Protection Set



High Protection Set



Medium Protection Set



Patient Protection Set



Basic Protection Set

Individually packed sets are configured to achieve the perfect level of isolation.

- The **operator and assistant set** includes the necessary disposables to wear during the entire day.
- The **patient set** protects the patient and also includes two gown covers for operator and assistant to wear for every new patient.



PREMIUM HIGH PROTECTION SET

FOR OPERATOR AND ASSISTANT*



Quantity	1 x	1 x	1 x	1 x
Description	Isolator™ N95 Particulate Respirator	SMS/PP Water resistant helmet-like cap with laces	Surgically folded Expo model gown (70 g) wrapped in medical paper & 2 hand towels	Face shield
Partcode	GPRN95	30.M0046.00	Operator: Size XL (length: 133 cm): 30.D1276.00 Assistant: Size S (length 118 cm): 30.D1278.00	10.M0200.00

Operator Set: **10.T0799.00** (non sterile) | **12.T0799.00** (sterile)
Assistant Set: **10.T0798.00** (non sterile) | **12.T0798.00** (sterile)



* The only difference between Operator Set and Assistant Set is the size of the gown.



HIGH PROTECTION SET

FOR OPERATOR AND ASSISTANT*



[MORE INFO](#)



[MORE INFO](#)



[MORE INFO](#)

Quantity	1 x	1 x	1 x
Description	Isolator™ N95 Particulate Respirator	SMS/PP Water resistant helmet-like cap with laces	Surgically folded Special Plus model gown (43 g) wrapped in medical paper & 2 hand towels
Partcode	GPRN95	30.M0046.00	Operator: Size XL (length: 118 cm): 30.D2173.00 Assistant: Size M (length: 112 cm): 30.D2171.00

Operator Set: **10.T0564.00** (non sterile) | **12.T0564.00** (sterile)

Assistant Set: **10.T0594.00** (non sterile) | **12.T0594.00** (sterile)

* The only difference between Operator Set and Assistant Set is the size of the gown.

**PACK
OF 10**



MEDIUM PROTECTION SET

FOR OPERATOR AND ASSISTANT*



[MORE INFO](#)



[MORE INFO](#)



[MORE INFO](#)

Quantity	1 x	1 x	1 x
Description	Ultra™ Sentitive Fogfree™ Earloop Mask with Secure Fit™ Mask Technology	Water resistant helmet-like cap with laces, SMS/PP, blue & white	Surgically folded Essential SMMS model gown (35 g)
Partcode	GCFCXSFSF	30.M0046.00	Operator: Size L (length: 125 cm): 30.D2302.00 Assistant: Size L (length: 125 cm): 30.D2302.00

Operator Set: **10.T0600.00** (non sterile) | **12.T0600.00** (sterile)

Assistant Set: **10.T0601.00** (non sterile) | **12.T0601.00** (sterile)

* The only difference between Operator Set and Assistant Set is the size of the gown.

**PACK
OF 10**



PATIENT PROTECTION SET FOR PATIENT

Each set has been designed to dress and protect the patient. It also contains 2 gown covers for operator and assistant to wear every new patient treatment.

<div> <div> <p>INCLUDES TWO GOWN COVERS TO PROTECT OPERATOR'S AND ASSISTANT'S GOWN.</p> <p>MORE INFO</p> </div> <div> <p>MORE INFO</p> </div> <div> <p>MORE INFO</p> </div> <div> <p>MORE INFO</p> </div> </div>				
Quantity	2 x	1 x	1 x	1 x
Description	Gown Cover - PE	Patient Safety Mask with shield	Round cap with elastic band	Sani-Tab™ Chain-Free Patient Towel
Partcode	30.D1002.00	GCPAT	30.M0040.00	WEXBLT
<div> <div> <p>MORE INFO</p> </div> <div> <p>MORE INFO</p> </div> <div> <p>MORE INFO</p> </div> <div> <p>MORE INFO</p> </div> </div>				
Quantity	2 x	1 x	1 x	2 x
Description	Patient Shoe Cover	Waterproof headrest with elastic band, round 24 cm Ø	Comfort Plus™ Premium Saliva Ejectors, White/Green Mint scented	2Cross20 Adhesive films 20 x 20 cm LDPE - BLUE
Partcode	30.D4003.00	20.O0011.00	ZWGCPM	32.U0020.00

Patient Set:

10.S9378.00 (non sterile)

12.S9378.00 (sterile)


PACK
OF 10


STAFF PROTECTION



BASIC PROTECTION SET FOR OPERATOR AND ASSISTANT*

The individually packed sets are composed in ideal configurations for short routine procedures.

		
Quantity	1 x	1 x
Description	Isolator™ N95 Particulate Respirator	Gown with elastic bands (length: 110 cm, 30 g), one size
Partcode	GPRN95	30.D1307.00

		
Quantity	1 x	1 x
Description	Isolator™ N95 Particulate Respirator	Shell Cover gown cover thumb up (length: 125 cm, 31 g) – one size only XL
Partcode	GPRN95	30.D1002.00

N95 & NWF Set:

10.T9118.00

PACK
OF 10

N95 & Shell Cover Set:

10.T9119.00

PACK
OF 10

* The only difference between Operator Set and Assistant Set is the type of gown.

PATIENT PROTECTION

OPERATOR PROTECTION

PRODUCT FEATURES OF EACH ITEM

Face Masks



Isolator™ N95 Particulate Respirator **GPRN95**

- NIOSH approved N95 Respirator.
- Meets guidelines for TB exposure control.
- Resistance to penetration by synthetic blood: 160 mmHg.
- Inner/outer fluid resistant layers.
- Chin pocket construction.
- EN 149: 2001 FFP2 NR / NIOSH N95

Although the N95 closely resembles a surgical mask, it is a respirator and must be used in accordance with all local regulations regarding respiratory protection.



Ultra™ Sentitive Fogfree™ Earloop Mask with Secure Fit™ Mask Technology **GCFCXSFSF**

- Type IIR Enhanced Fluid Resistance.
- ASTM 3.
- Fluid resistant outer layer.
- Fluid resistant, extra soft white hypoallergenic inner cellulose layer – will not lint, tear or shred.
- Free of chemicals, inks and dyes.
- Extra-long aluminum nose piece.
- FogFree™ Strip on mask inner layer.
- Secure Fit™ Mask Technology



Patient Safety Mask with shield **GCPAT**

- Fluid resistant. Covers patient's eyes and nose.
- Protects fragile eye and mucosal tissue from flying particles, spray and splatter.
- Makes patients feel safe, but not "in the dark," with attached clear-vision shield.
- Form-fitting profile, flexible materials provide maximum access to mouth.
- White spunbond inner and outer layers.



Face shield **10.M0200.00**

- Wide polyethylene fog preventing shield.
- Polyurethane insert.
- Suitable for people wearing glasses.

Caps



SMS/PP Water resistant helmet-like cap with laces **30.M0046.00**

- Water / Fluid resistant - Breathable part.



Round cap with elastic band (for the patient) **30.M0040.00**

- Enables ventilation.
- Regarding respiratory protection.

Gowns



Expo model gown (Sontara®, 70 g), surgically folded, wrapped in medical paper with 2 hand towels

30.D1276.00 (XL, length: 133 cm) /
30.D1278.00 (S, length: 118 cm)

- Spunlace medical nonwoven fabric.
- High liquid-repellence.
- Extraordinary sense of comfort.
- Velcro closure, cotton cuffs.
- Double internal and external side laces.
- Surgically folded.
- 2 hand towels.
- Wrapped in medical paper.



Surgically folded SMMS model gown (35 g) **30.D2302.00** (L, length: 125 cm)

- Water-repellent and breathable, SMMS.
- Cotton cuffs.
- Waist-belt.
- Surgical folding make.
- Ultrasonic welded seams to guarantee a more efficient barrier against bacteria and fluids.



Special Plus model gown (SMS Microfibre, 43 g), surgically folded, wrapped in medical paper with 2 hand towels

30.D2173.00 (XL, length: 118 cm) /
30.D2171.00 (M, length: 112 cm)

Water-repellent and breathable, SMS Microfibre is used to produce these gowns providing an excellent relationship between comfort and performance. Its cotton cuffs, side-fastened waist-belt, surgical folding and internal towel make it the ideal choice for clinical work sessions.



Gown with elastic bands (30 g), (NWF)
30.D1307.00 (one size, length: 110 cm)

With elastic bands, ideal for short routine procedures or as a supplement to patient protection and isolation from pathogens.



Shell Cover gown cover thumb up (31 g)
30.D1002.00 (one size only XL, length: 125 cm)

Thumb Up, perfect for routine practices: the special Medical-grade PE guarantees a complete impermeability to pathogenic and contaminant agents, while protecting the professional uniform from any accidental splash or contamination. The additional Thumb-Up system helps keeping the gown fixed while worn.

- Medical grade
- Waterproof
- High insulation standard

Miscellaneous



Sani-Tab™ Chain-Free Patient Towel
WEXBLT

- Easy to use – simply position towel anywhere on the patient.
- Bib stays where placed to maximize coverage, unlike bib clips and other .
- Disposable towel holders.
- Can be placed high on patient to protect clothing.



Patient Shoe Cover
30.D4003.00

- Waterproof shoe s with elastic opening.



Waterproof headrest with elastic band – white, round 24 cm Ø
20.00011.00

- Waterproof and absorbent NWF protections.



Comfort Plus™ Premium Saliva Ejectors
White/Green Mint scented
ZWGCPM

- The kinder, gentler saliva ejector...takes in fluid, not tissue!
- Patented tip design enhances patient comfort – will not grab mucosal tissue.
- Facilitates uninterrupted fluid and light debris removal.
- Optimizes fluid removal as patented tip design prevents tissue from blocking suction.
- Holds shape when formed as result of proprietary manufacturing process.
- Respiratory protection.



2Cross20 Adhesive films
20 x 20 cm LDPE - BLUE
32.U0020.00

- LDPE protections.
- Avoid cross-contamination.



STERILE

OMNIA DRESSING SETS

Simplify your dental routine with sterile Dressing Sets.



OPERATOR SET

This sterile set has been designed to dress the operator. It is ideal for all dental implant surgery, parodontology and maxillofacial operations. The set is supplied in a double sterile packaging, wrapped in medical paper. The practical peel-off label allows easy traceability of the operation.

Quantity	1 x	1 x	1 x	2 x
Description	Surgically folded Expo model gown (size L) and 2 hand towels wrapped in medical paper	Surgical cap with adjustable strings	SPIRITUS PLUS® face mask with strings	Non-slip shoe covers
Material	Sontara®	NWF	NWF	NWF
Color	Light blue	Light blue	Light blue	Light blue

Operator Set
12.T0010.00 (sterile)



“ The ability to provide bespoke drape and gown packs for both Specialist surgeons means that we can continue to maintain the quality of care necessary for our patients. Omnia were very helpful with setting up the different packs and their service has been excellent throughout the whole process.

Dr. Paul Stone BDS, FDS, RCS(Edin)
Specialist Oral Surgeon & Clinical Director, Blackhills Clinic - UK



ASSISTANT SET

This sterile set has been designed to dress the assistant and the visiting staff in the operating theatre. It is ideal for all dental implant surgery, parodontology and maxillofacial operations. The set is supplied in a double sterile packaging, wrapped in medical paper, and the gown is surgically folded. The practical peel-off label allows easy traceability of the operation.

				
Quantity	1 x	1 x	1 x	2 x
Description	Surgically folded Special Plus model gown (size M) and 2 hand towels wrapped in medical paper	Unisex breathable cap	SPIRITUS PLUS® face mask with strings	Non-slip shoe covers
Material	SMS Microfibre	NWF	NWF	NWF
Color	Light blue	Light blue	Light blue	Light blue

Assistant Set

12.T0001.00 (sterile)



PATIENT SET

This sterile set has been designed to dress the patient. It is ideal for all dental implant surgery, parodontology and maxillofacial operations. The set is supplied in a double sterile packaging, wrapped in medical paper. The gown has elastic bands at the wrists for a better comfort of the patient. The practical peel-off label allows easy traceability of the operation.

			
Quantity	1 x	1 x	1 x
Description	Fluid-repellent gown with elastic	Round cap with elastic band	Liquid proof shoe covers
Material	NWF	NWF	NWF
Color	White	White	Light blue

Patient Set

12.T0002.00 (sterile)



OPERATOR KIT

This sterile kit has been designed to dress the operator. It is ideal for all not very complicated or short duration operations. The practical peel-off label allows easy traceability.

				
Quantity	1 x	1 x	1 x	2 x
Description	Fluid-repellent gown with cotton cuffs	Surgical cap with adjustable strings	SPIRITUS PLUS® face mask with elasticated ear loops	Liquid proof shoe covers
Material	NWF	NWF	NWF	PVC
Color	Light blue	Light blue	Light blue	Light blue

Operator Kit

12.T0107.00 (sterile)



ASSISTANT KIT

This sterile kit has been designed to dress the assistant or visiting staff. It is ideal for all not very complicated or short operations. The practical peel-off label allows easy traceability of the operation.

				
Quantity	1 x	1 x	1 x	2 x
Description	Fluid-repellent gown with cotton cuffs	Round cap with elastic band	SPIRITUS PLUS® face mask with elasticated ear loops	Liquid proof shoe covers
Material	NWF	NWF	NWF	PVC
Color	Light blue	Light blue	Light blue	Light blue

Assistant Kit

12.T0108.00 (sterile)



ORAL SURGERY SET

Sterile set of drapes, ideal for all minor operations, such as wisdom teeth extractions and second-stage surgeries. The special self-adhesive U-shaped cut on one side of the drape allows circumscribing the surgical site effectively. The practical peel-off label allows easy traceability.



Quantity	1 x	2 x
Description	Absorbent/liquid proof drape 50 x 75 with adhesive U-shaped cut 11 x 9 cm	Absorbent/liquid proof drapes 50 x 75 cm
Material	NWF + PE	NWF + PE
Color	Light blue	Light blue

Oral Surgery Set
12.T0105.00 (sterile)

**PACK
OF 10**



I like Omnia products because of quality and reliability. I am using the individual kits for all my implant and bone regeneration surgery. The products make our surgical set-up infection resistant.

Istvan Urban DMD, MD, PhD
DMD, MD, PhD Associate Professor University of Szeged,
Department of Periodontology Szeged, Hungary



I remember when we used to use long fabric drapes to cover the patient. Luckily it is an old memory.

Thanks to Omnia products, we are more efficient in time management and always up to standards.

Mario Kirste
MSc Implantologie Frankfurt/Oder



BASIC IMPLANT SET

The Basic Implant Set has been designed for all routine surgical activities where a surgical protocol ensuring an adequate degree of sterility and guaranteeing the safety of operators and patients is needed. Due to its reasonable price, it is especially fit for all surgical procedures with immediate loading or for operations, such as re-opening of implants, flap surgery, operations on minor surgical wounds.



Quantity	1 x	1 x	2 x	1 x
Description	Surgically folded Special Plus model gown (size L) and 2 hand towels wrapped in medical paper	Gown with cotton cuffs	Surgical caps with fixed laces	Round surgical cap for the patient
Material	SMS Microfibre	NWF	Absorbent NWF	Absorbent NWF
Color	Light blue	Light blue	Light blue	Light blue



Quantity	1 x	1 x	1 x	1 x
Description	Absorbent/liquid proof drape 50 x 75 cm with self-adhesive* side	Absorbent/liquid proof drape 75 x 90 cm with adhesive* U-shaped cut 11x9 cm	Omnisleeve sheath 120 x 7 cm	Disposal bag for contaminated material
Material	NWF + PE	NWF + PE	PE	PVC
Color	Light blue	Light blue	Transparent	Red

Basic Implant Set

12.S5621.00 (sterile)

* Bold numbers indicate the adhesive side



IMPLANTOLOGY SET

Complete sterile set including all items needed to perform any kind of oral implantology operation. It contains the gowns for both the operator and the assistant, surface drapes with adhesive fixing strips, sheaths for hand pieces, films for handles and keypads, patient drape with self-adhesive U-shaped. The set is wrapped in an absorbent/liquid proof drape 145 x 160 cm.



Quantity	2 x	2 x	1 x	1 x
Description	Surgically folded Special Plus model gowns (size L) and hand towel wrapped in medical paper	Liquid proof/absorbent drapes with self-adhesive* side 75 x 90 cm	Surgical aspirator with Yankauer tip	Liquid proof/absorbent drape 100x150 cm with adhesive* U-shaped cut 6,5 x 30 cm
Material	SMS Microfibre	NWF + PE	PVC	NWF + PE
Color	Light blue	Light blue	Transparent	Light blue



Quantity	1 x	2 x	1 x
Description	Omnislave sheath 120 x 7 cm	2Cross20 adhesive films 20 x 20 cm	Disposal bag for contaminated material
Material	PE	LDPE	PVC
Color	Transparent	Light blue	Red




Implantology Set

12.S4671.00 (sterile)



ORAL SURGERY SET DRAPE WITH ADHESIVE U-SHAPED

Complete sterile set for complex procedures, such as important rehabilitations, bone distractions, maxillofacial operations. It contains all what is needed to prepare the operating theatre for demanding operations: surface drapes with adhesive fixing strips, **breathable drape for the patient made of Sontara® with self-adhesive U-shaped cut**, sheaths for hand pieces, adhesive films for keypads, special cover in medical PVC for Physiodispenser, set for patient disinfection, Mayo table cover, 2 surgical gowns in Sontara® with 2 sterile hand towels. The set is wrapped in a practical absorbent, liquid rproof drape 160x190 cm and packed in a Tyvek® envelope.

					
Quantity	2 x	1 x	2 x	1 x	1 x
Description	Surgically folded Expo model gown (size L) and 2 hand towels wrapped in medical paper	Round cap with elastic band	Liquid proof/ absorbent drapes 75 x 90 cm with adhesive* side	Mayo table cover 80 x 120 cm	Cover 40 x 50 cm with adhesive* side
Material	Sontara®	NWF	NWF + PE	NWF + PVC	PE
Color	Light blue	Light blue	Light blue	Light blue	Transparent
					
Quantity	2 x	1 x	2 x	1 x	1 x
Description	Omnisleeve sheaths 240 x 7 cm	ASP51/F surgical aspirator	2Cross20 adhesive films 20 x 20 cm	Sponge with grip for patient disinfection	Bowl 250 ml
Material	PE	PVC	LDPE	PVC	PE
Color	Transparent	Transparent	Light blue	Blue	Blue
					
Quantity	1 x	1 x	5 x	1 x	
Description	Fluid-repellent drape 67 x 75 cm with adhesive* side	Fluid-repellent drape 133x200 cm with adhesive* U-shaped 6,5 x 50 cm and loop	Round tampons Ø 3 cm	Disposal bag for contaminated material	
Material	Sontara®	Sontara®	Cotton	PVC	
Color	Light blue	Light blue	White	Red	

ORAL SURGERY SET – DRAPE WITH ADHESIVE U-SHAPED

12.S4041.00 (sterile)














* Bold numbers indicate the adhesive side



ORAL SURGERY SET ADHESIVE HOLE SURGICAL DRAPE**

** for those who prefer working, using a surgical protocol that envisages draping the patient's eyes

Complete sterile set for complex procedures, such as important rehabilitations, bone distractions, maxillofacial operations. It contains all what is needed to prepare the operating theatre: surface drapes with adhesive fixing strips, **breathable, self-adhesive fenestrated drape for the patient made of Sontara®**, sheaths for hand pieces, adhesive films for keypads, special cover in medical PVC for Physiodispenser, set for patient disinfection, Mayo table cover, 2 surgical gowns in Sontara® with 2 sterile hand towels. The set is wrapped in a practical absorbent, liquid proof drape 160x190 cm and packed in a Tyvek® envelope.

					
Quantity	2 x	1 x	2 x	1 x	1 x
Description	Surgically folded gowns (size L), 2 hand towels wrapped in medical paper	Round cap with elastic band	Liquid proof/absorbent drapes 75 x 90 cm with self-adhesive* side	Mayo table cover 80 x 120 cm	Cover 40 x 50 cm with adhesive* side
Material	Sontara®	NWF	NWF + PE	NWF + PVC	PE
Color	Light blue	Light blue	Light blue	Light blue	Transparent
					
Quantity	2 x	1 x	2 x	1 x	1 x
Description	Omnisleeve sheaths 240 x 7 cm	ASP51/F surgical aspirator	2Cross20 adhesive films 20 x 20 cm	Sponge with grip for patient disinfection	Bowl 250 ml
Material	PE	PVC	LDPE	PVC	PE
Color	Transparent	Transparent	Light blue	Blue	Blue
					
Quantity	1 x	5 x		1 x	
Description	Fluid-repellent drape 133 x 200 cm with adhesive* off-centre hole 6 x 9 cm and loop	Round tampons Ø 3 cm		Disposal bag for contaminated material	
Material	Sontara®	Cotton		PVC	
Color	Light blue	White		Red	

Oral Surgery Set – Adhesive Hole Surgical Drape

12.S4096.00 (sterile)


* Bold numbers indicate the adhesive side




ORAL SURGERY SET NO PANIC DRAPE

Complete sterile set for complex procedures, such as important rehabilitations, bone distractions, maxillofacial operations. It contains all what is needed to prepare the operating theatre: surface drapes with adhesive fixing strips, breathable drape for the patient made of Sontara® with a **protection window made of a transparent material that allows an easy and immediate communication between the patient and the operator, avoiding unpleasant panic situations**, sheaths for hand pieces, adhesive films for keypads, special cover in medical PVC for Physiodispenser, set for patient disinfection, Mayo table cover, 2 surgical gowns in Sontara® with 2 sterile hand towels each.

The set is wrapped in a practical absorbent, liquid proof drape 160 x 190 cm and packed in a Tyvek® envelope.



Quantity	2 x	1 x	1 x	2 x
Description	Expo model gowns (size L), surgically folded with 2 hand-towels wrapped in medical paper	Liquid proof gown with cuffs	Round cap with elastic band	Absorbent/water-proof drapes 75 x 90 cm with adhesive* side
Material	Sontara®	Non woven	Non woven	NWF+ PE
Color	Blue	Blue	Blue	Blue



Quantity	1 x	1 x	2 x	1 x
Description	Mayo's drape 80 x 120 cm	PVC Transparent drape 40 x 50 cm with adhesive* side	Tubing sleeves 240 x 7 cm with Omnisleeve inserter	Surgical aspirator ASP51/F
Material	NWF+PVC	PVC	PE	PVC
Color	Blue	Transparent	Transparent	Transparent



Quantity	2 x	1 x	1 x	1 x
Description	Adhesive films 21 x 31 cm	Manual irrigation syringe 10 cc	Sponge for patient disinfection	Bowl 250 ml
Material	PVC	PP	PE+PU	PE
Color	Blue	Transparent	Blue	Blue



Quantity	1 x	5 x	4 x	1 x
Description	Liquid repellent surgical drape cm 133 x 200 with triangular shape adhesive* 11 x 10 cm in the low side, trans- parent opening for the eyes and velcrum loop to fix tubing	Gauzes sponges Ø 3 cm	Folded compresses 5 x 5 cm, 4-ply	Waste bag for contaminated material
Material	Sontara®	Cotton	Non woven	PVC
Color	Blue	White	White	Red

Oral Surgery Set – No panic drape
12.S4061.00 (sterile)

* Bold numbers indicate the adhesive side



PARODONTALSET

Sterile set for periodontal surgery or small oral surgery. It contains patient drapes and surface drapes. It includes sheaths for hand pieces, adhesive films for handles and keypads and folded compresses.

The set is packed in a double packaging. The practical peel-off label allows easy traceability of the operation.

Quantity	1 x	1 x	2 x	2 x	2 x
Description	Adhesive* U shaped 11 x 9 cm absorbent/water-proof drape 75 x 90 cm	Absorbent/water proof drape 75 x 90 cm	Absorbent/water proof drapes 50x75 cm	Tubing sleeves 120 x 7 cm with fixing adhesives	Folded compresses 5 x 5 cm, 4-ply
Material	Non woven + PE	Non woven + PE	Non woven + PE	LDPE	Non woven
Color	Blue	Blue	Blue	Transparent	White

Parodontalset

12.T0013.00 (sterile)



DENTSET

Sterile set containing patient drapes and surface drapes. It includes sheaths for hand pieces and adhesive films for handles and keypads. The set is packed in a double packaging. The practical peel-off label allows easy traceability of the operation.

Quantity	1 x	1 x	2 x	2 x	4 x
Description	Absorbent/liquid proof drape 100 x 150 cm with U-shaped cut 6,5 x 30 cm*	Absorbent/liquid proof drape 75 x 90 cm	Absorbent/liquid proof drapes 50 x 75 cm	Omnisleeve sheaths 120 x 7 cm	2Cross20 adhesive films 20 x 20 cm
Material	NWF + PE	NWF + PE	NWF + PE	PE	LDPE
Color	Light blue	Light blue	Light blue	Transparent	Light blue

Dentset

12.T0003.00 (sterile)

* Bold numbers indicate the adhesive side



“ Omnia products has given me an efficient, safe and high quality way to control sterility for my patients during surgical implant procedures.

Dr. Sascha A. Jovanovic
Periodontist and Implant Surgeon, Los Angeles - California



“ The packs for the operating room and patient preparation allow optimal control of sterility, an important cornerstone (...). In addition, they are extremely practical and easy to use for the whole team.”

Prof. Tiziano Testori MD DDS FICD,
Doctor-Surgeon specialist in dentistry
and orthodontics – Como, Italy

CUSTOMIZED PROCEDURAL SETS

To satisfy your needs.

The painstaking, rigorous and proper preparation of the surgical environment following a pre-arranged protocol is an essential element for a successful outcome of the operation. The important task consists of ensuring the availability, choosing, checking and preparing the large number of individual products needed for a surgical procedure.

This process is extremely time-consuming for operating theatre professionals. To satisfy the needs of the operating theatre staff, OMNIA® developed a solution that is extremely efficient: the **customized Procedural Sets**.

A single Procedural Set containing all the individual items that are used for preparing **your procedures in your operating theatre**.

CUSTOMIZATION

A customized Procedural Set contains all disposable medical devices prepared according to your particular requirements, in order to meet the needs of your individual surgical protocols. The aim is to improve the efficiency of the operating theatre, by simplifying the preparation, reducing preparation delays, enhancing performance quality and the quality of assistance to patients.

With optimized time and cost it allows you to focus more on the clinical case.

PRODUCTS FOR YOUR CUSTOMIZED PROCEDURAL SET

We offer a wide range of products with many variations to suit your specific needs:

- PPE
- Surgical gowns
- Patient and surface drapes
- Surgical accessories (sheaths and adhesive films)
- Aspiration systems
- Irrigation systems
- and more...

On our website you will find a range of products which are not present in the catalogue; these can only be used to create a personalized set.

TIPS FOR CREATING YOUR CUSTOMIZED PROCEDURAL SET

- List the items you would like in your set
- Choose the materials, quantities and measurements of the items
- Estimate the number of procedures planned during the year
- Sample and (if necessary) modify your choice of products.



CLICK HERE
TO VIEW OUR WEBSITE!

THE EASY WAY TO YOUR CUSTOMIZED PROCEDURAL SET

MODULE

Tell us what products you want in the set: sizes, quantities and the exact order you want to find them.



SAMPLE SET

Our product specialists assemble your Customized Procedural Set Sample.



SAMPLING

Receive your sample set & test in your practice, with your staff if it meets all your requirements.



PRODUCTION

After your first order follows the productions of your Customized Procedural Set.**



DELIVERY

After approx. 4 weeks, you will receive your Customized Procedural Sets in your practice.



OUR SERVICE: **PERSONAL CONSULTATION**

Our network of consultants is available to guide you so you can achieve the best Customized Procedural Set for your needs. You can contact us directly, or our highly trained agents and dealer partners. Together, we help you to choose the ideal components and features in order to create the perfect Customized Procedural Set for you.

OUR SERVICE: **COMPLIANCE AND TRACEABILITY**

Complies with all required standards* that regulate production specifications and the safe placing on the market of the devices, as well as their preservation, transport, labelling, after-sale surveillance, management of claims and non-conformities and, in particular, the traceability of the devices.

The batch number assigned to the device, its code and expiry date, are the elements identifying the product and guaranteeing its traceability. The labels used by Omnia contain all the information needed for the identification of the product. It has two detachable parts for the end user allowing an easy data storage and traceability.

*Medical Device Directive 93/42/EEC (updated by 2007/47/EC) and UNI EN ISO 13485:2004 standards. **Minimum order quantity: 30 pcs.



For more information on Hu-Friedy products or to place an order, call 00800 48 37 43 39 or visit Hu-FriedyGroup.eu

For more information on Crosstex products, call +39 05 24 52 74 53 or visit Crosstex.com

For more information on Omnia products or to place an order, call +39 05 24 52 74 53 or visit OmniaSpa.eu

Hu-Friedy Mfg. Co., LLC • European Headquarters • Lyoner Str. 9 • 60528 Frankfurt am Main, Germany • HuFriedyGroup.eu

As some of the products included in this catalogue might not be available in some countries, we invite you to verify the availability of the product with your local authorised dealer or to contact our customer care at the free phone number 00800 48 37 43 39. Sontara® is a trademark of Sontara® South Asia. and Tyvek® is trademark of DuPont Safety & Construction, Inc.. All product names are trademarks of Omnia Srl, a Cantel Medical Company, its affiliates, or related companies, unless otherwise noted

©2021 Omnia Srl. All rights reserved. 90.L0093.00/1121