

# FACE MASKS AND FACE PROTECTION



# CROSSTEX MASKS

**BETTER FIT.  
BETTER FILTRATION.  
BETTER PROTECTION.**



**CHIN CLOSURE**

## ADDITIONAL AVAILABLE OPTIONS:

- **FOG-FREE™** FogFree™ Strip blocks and absorbs moisture; forms a strong seal preventing fogging of eyewear while it cushions the nose ridge.
- **ULTRA SENSITIVE** Fluid resistant, extra soft white hypoallergenic inner cellulose layer.
- **BIOSAFE®** Features a long-lasting protective bactericidal coating that is permanently bound to the mask and is effective against many harmful bacteria.

**Award-winning face masks.** 

The first to feature adjustable aluminum nose and chin closures, providing a safe, secure fit, regardless of face shape or size.

## **SECUREFIT™ MASK TECHNOLOGY**

Proven to provide **better protection** than standard earloop masks.<sup>1</sup>

Snug, custom fit minimizes dangerous gaps for superior protection.

# PERSONAL PROTECTION EQUIPMENT (PPE)

## THE HEART OF ANY INFECTION PREVENTION PROGRAM.

Crosstex manufactures an extensive portfolio of PPE to meet all your practice needs, including our Made-in-USA and Made in Italy face masks that meet EN 14683 Type and ASTM Level classifications. Selecting the appropriate mask for a particular procedure is a critical component of your Personal Protective Equipment (PPE) protocol. Although masks may look similar, each mask has notable differences affecting its level of protection. Understanding EN14683 and ASTM performance can help ensure your mask will provide appropriate protection to minimize the spread of potentially infectious diseases.



EN14683 Classification	TYPE IIR (Enhanced Fluid Resistance, ASTM Level 3)	TYPE IIR	TYPE II
Procedures	Ideal for procedures where heavy to moderate amounts of fluid, spray and/or aerosols are produced. <ul style="list-style-type: none"> <li>• Complex oral surgery</li> <li>• Crown preparation</li> <li>• Implant placement</li> <li>• Periodontal surgery</li> <li>• Use of ultrasonic scalers</li> <li>• Laser based applications<sup>2</sup></li> </ul>	Ideal for procedures where moderate to light amounts of fluid, spray and/ or aerosols are produced. <ul style="list-style-type: none"> <li>• Limited oral surgery</li> <li>• Endodontics</li> <li>• Prophylaxis</li> <li>• Restorative</li> <li>• Sealant</li> </ul>	Ideal for procedures where light amounts of fluid, spray and/or aerosols are produced <ul style="list-style-type: none"> <li>• Patient exams</li> <li>• Operatory cleaning</li> <li>• Impressions</li> <li>• Lab finishing &amp; polishing</li> <li>• Orthodontics</li> </ul>
Fluid resistance, kPa (mm Hg)	21.3 (160)	16.0 (120)	Not Required
BFE, @ 3 microns	≥98%	≥98%	≥98%
Delta P, Pa/cm <sup>2</sup> (mm H <sub>2</sub> O/cm <sup>2</sup> )	< 49.0 (<5)	< 49.0 (<5)	< 29.4 (<3)
Microbial Cleanliness, cfu/g	< 30	< 30	< 30
<b>Additional ASTM 2100 Classification</b>	<b>Level 3</b>	<b>Level 2</b>	<b>Level 1</b>
PFE, @ 0.1 micron	≥98%	≥98%	≥95%
Flame spread	Class 1	Class 1	Class 1
Fluid resistance, kPa (mm Hg)	21.3 (160)	16.0 (120)	10.7 (80)
Delta P, Pa/cm <sup>2</sup> (mm H <sub>2</sub> O/cm <sup>2</sup> )	< 49.0 (<5)	< 39.2 (<4)	< 29.4 (<3)
Product families meeting both EN14683 & ASTM classifications	Ultra™ Masks Ultra™ Sensitive Masks Ultra™ Sensitive Masks with BIOSAFE®	Procedural Surgical Tie-On	Isofluid™ Masks Isofluid™ Plus Masks Le Petit Mask

### HOW OFTEN SHOULD I CHANGE MY MASK?<sup>2</sup>



Non-aerosol  
generating procedures  
**every 60 minutes**



Aerosol  
generating procedures  
**every 20 minutes**



**Every each  
patient**



If the mask gets  
**soiled or wet**



**FILTRATION LEVEL KEY** ● EN 149: 2001 FFP2 NR ● Type IIR Enhanced Fluid Resistance ASTM 3 ● Type IIR ASTM 2 ● Type II ASTM 1 ● Minimum Performance

<sup>2</sup> Masks are considered a secondary control and are not meant to replace recommended primary engineering controls for laser plume exposure.

CAUTION: Like patients, every procedure is unique. Use your clinical judgment to determine the appropriate level of barrier protection based on the length of the procedure, the amount of fluid or aerosol, and standard precautions.



## MASKS STANDARDS

### EN14683 & ASTM 2100 specification definitions

- **Fluid resistance** measures the mask's ability to prevent fluids from traveling through the material. The higher the fluid resistance (filtration), the better the protection. (Measured in kPa or mm Hg.)
- **BFE (bacterial filtration efficiency)** measures the mask's ability to prevent bacteria from traveling through the material. (Measured at 3 microns.)
- **Delta P (differential pressure)** determines breathing resistance—the higher the Delta P, the less the breathability, but the better the filtration. (Measured in Pa/cm<sup>2</sup> or mm H<sub>2</sub>O/cm<sup>2</sup>.)

### Additional EN14683 specification definition

- **Microbial Cleanliness** measures the total bioburden of the mask per gram of mask material tested. (Measured in colony forming units per gram; cfu/g.)


### Additional ASTM 2100 specification definitions

- **PFE (particulate filtration efficiency)** measures the mask's ability to prevent sub-micron particulates from traveling through the material. The size of the particles filtered is critical. (Measured at 0.1 micron.)
- **Flame spread** is a ranking derived by laboratory standard test methodology of a material's propensity to burn rapidly and spread flames.

*Masks are considered a secondary control and are not meant to replace recommended primary engineering controls for laser plume exposure.*

*CAUTION: Like patients, every procedure is unique. Use your clinical judgment to determine the appropriate level of barrier protection based on the length of the procedure, the amount of fluid or aerosol, and standard precautions.*

# MASK FEATURE GUIDE



- A** Extra-long, enclosed aluminum nose closure forms a strong seal around nose and cheekbones for maximum protection.
- B** Optional attached wrap-around face shield is optically clear and distortion-free.
- C** The elastic, flat ear loops minimize pressure to the ears and are attached to the outside of the mask to eliminate skin irritation and provide the most comfortable fit.
- D** Enclosed aluminum chin closure significantly reduces gapping at the sides and bottom of the mask, minimizing the potential of inhaling airborne contaminants.
- E** Unique vapor barrier on the outside layer blocks moisture and prevents fogging.
- F** Fluid resistant, spunbond outer layer for maximum protection.
- G** Ultrasoft, odorless, medical grade facial tissue outer layer.
- H** Highest quality meltblown filter media for best filtration efficiency and breathability.
- I** Extra soft, white hypoallergenic inner cellulose layer for sensitive skin is free of chemicals, inks and dyes; will not lint tear or shred for ultimate comfort.
- J** Fluid resistant, white, spunbond inner layer.
- K** Soft, white medical-grade facial tissue inner layer.
- L** FogFree™ Strip blocks and absorbs moisture; forms a strong seal preventing fogging of eyewear while it cushions the nose ridge.
- M** Sof-Bond™ Ultrasonic Mask assembly ensures strongest construction and highest quality, eliminating holes and defects.

	A	B	C	D	E	F	G	H	I	J	K	L	M
Ultra™ Sensitive Earloop Mask with Secure Fit™ Mask Technology	•		•	•		•		•	•				•
Ultra™ Sensitive FogFree™ Earloop Mask with Secure Fit™ Mask Technology	•		•	•	•	•		•	•			•	•
Ultra™ Sensitive FogFree™ Earloop Mask with shield with Secure Fit™ Mask Technology	•	•	•	•	•	•		•	•			•	•
Ultra™ Earloop Mask with Secure Fit™ Mask Technology	•		•	•		•		•		•			•
Ultra™ Sensitive Earloop Mask	•		•			•		•	•				•
Ultra™ Sensitive FogFree™ Earloop Mask	•		•		•	•		•	•			•	•
Ultra™ Sensitive FogFree™ Earloop Mask with shield	•	•	•		•	•		•	•			•	•
Ultra™ Sensitive Earloop Mask treated with BIOSAFE™	•		•			•		•		•			•
Ultra™ Earloop Mask	•		•			•		•		•			•
Ultra™ FogFree™ Earloop Mask	•		•		•	•		•		•		•	•
Ultra™ FogFree™ Earloop Mask with shield	•	•	•		•	•				•		•	•
Procedural Earloop Mask with Secure Fit™ Mask Technology	•		•	•		•		•		•			•
Procedural Earloop Mask	•		•			•		•					•
Surgical Tie-on Mask	•		•			•		•					•
Isofluid™ Plus Pure Earloop Mask with Secure Fit™ Mask Technology	•		•	•		•		•	•				•
Isofluid™ Earloop Mask with Secure Fit™ Mask Technology	•		•	•		•		•			•		•
Isofluid™ FogFree™ Earloop Mask with Secure Fit™ Mask Technology	•		•	•		•		•			•	•	•
Isofluid™ FogFree™ Earloop Mask with shield with Secure Fit™ Mask Technology	•	•	•	•		•		•			•	•	•
Isofluid™ Plus Earloop Mask	•		•					•	•				•
Isofluid™ Earloop Mask	•		•					•		•			•
Isofluid™ FogFree™ Earloop Mask	•		•					•		•		•	•
Isofluid™ FogFree™ Earloop Mask with Shield	•	•	•					•		•		•	•
Isolite™ Earloop Mask	•		•				•	•			•		•

All Crosstex face masks are not made with natural rubber latex and are not reusable.

MASK ORDERING INFORMATION

PRODUCT DESCRIPTION, QUANTITY	QTY	CODE, COLOR		
ISOLATOR™ MASKS				
EN 149: 2001 FFP2 NR / NIOSH N95				
Isolator™ Plus N95 Particulate Respirator	28/Box	GPRN95, Blue/White Stripe		
ULTRA™ SENSITIVE, ULTRA™, BIOSAFE® MASKS				
EN 14683 Type IIR (enhanced fluid resistance) / ASTM Level 3				
Ultra™ Sensitive Earloop Mask with Secure Fit™ Mask Technology	50/Box	GCFCXSSF, White		
Ultra™ Sensitive FogFree™ Earloop Mask with Secure Fit™ Mask Technology	40/Box	GCFCXSFSF, White		
Ultra™ Sensitive FogFree™ Earloop Mask with shield with Secure Fit™ Mask Technology	25/Box	GCPWSSF, White		
Ultra™ Earloop Mask with Secure Fit™ Mask Technology	50/Box	GCFCXUSF, Blue		
Ultra™ Sensitive Earloop Mask	50/Box	GCFCXS, White		
Ultra™ Sensitive FogFree™ Earloop Mask	40/Box	GCFCXSF, White		
Ultra™ Sensitive FogFree™ Earloop Mask with shield	25/Box	GCPWS, White		
Ultra™ Sensitive Earloop Mask with BIOSAFE®	40/Box	GCFCXSBS, White		
Ultra™ Earloop Mask	40/Box	GCFCXU, Blue		
Ultra™ FogFree™ Earloop Mask	40/Box	GCFCX, Blue		
Ultra™ FogFree™ Earloop Mask with shield	25/Box	GCPW, Blue		
PROCEDURAL MASKS				
EN 14683 Type IIR / ASTM Level 2				
Procedural Earloop Mask with Secure Fit™ Mask Technology	50/Box	GCPBLSFEU, Blue	GCPLVSF, Lavender	GCPPKSF, Pink
Procedural Earloop Mask	50/Box	GCPBLEU, Blue GCPLV, Lavender	GCPPK, Pink GCPYE, Yellow	
Procedural Earloop Mask	100/Box	GCPBL100, Blue		
Surgical Tie-On Mask	50/Box	GCS, Blue		
ISOFLUID™ MASKS				
EN 14683 Type II / ASTM Level 1				
Isofluid™ Plus Pure Mask with Secure Fit™ Mask Technology	50/Box	GPLUSWHSF, White		
Isofluid™ Earloop Mask with Secure Fit™ Mask Technology		GCIBLSF, Blue	GCILVSF, Lavender	GCIPKSF, Pink
Isofluid™ FogFree™ Earloop Mask with Secure Fit™ Mask Technology	40/Box	GCICXBSF, Blue		
Isofluid™ FogFree™ Earloop Mask with shield with Secure Fit™ Mask Technology	25/Box	GCIPWBSF, Blue		
Isofluid™ Plus Earloop Mask	50/Box	GPLUSBL, Blue	GPLUSKA, Kaleidoscope	GPLUSPK, Pink
Isofluid™ Earloop Mask	50/Box	GCIBL, Blue GCIGR, Green GCILV, Lavender	GCIPK, Pink GCISA, Sapphire GCITE, Teal	GCITQ, Turquoise
Isofluid™ Earloop Mask	100/Box	GCIBL100, Blue		
Isofluid™ FogFree™ Earloop Mask	40/Box	GCICXB, Blue	GCICXS, Sapphire	GCICXT, Turquoise
Isofluid™ FogFree™ Earloop Mask with shield	25/Box	GCIPWB, Blue	GCIPWS, Sapphire	GCIPWT, Turquoise
Le Petit Masque	50/Box	GCAPTPK, Pink		
ISOLITE™ MASKS				
Minimum Performance				
Isolite™ Earloop Mask	50/Box	GCLBL, Blue	GCLPK, Pink	

# ULTRA™ SENSITIVE, ULTRA™, BIOSAFE® MASKS CROSSTEX

## TYPE IIR (Enhanced Fluid Resistance)

Ideal for procedures where heavy to moderate amounts of fluid, spray and/or aerosols are produced.

### ULTRA™ SENSITIVE EARLOOP MASK with Secure Fit™ Mask Technology

- Fluid resistant outer layer.
- Fluid resistant, extra soft white hypoallergenic inner cellulose layer – will not lint, tear or shred.
- Free of chemicals, inks and dyes.
- Extra-long aluminum nose piece.



#### ULTRA™ SENSITIVE EARLOOP MASK with Secure Fit™ Mask Technology

CODE	DETAILS	QTY
<b>GCFCXSSF</b>	White	50/Box

### ULTRA™ SENSITIVE FOGFREE™ EARLOOP MASK with Secure Fit™ Mask Technology

- Provides all the features of Ultra™ Sensitive Earloop Mask with Secure Fit™ Mask Technology with the addition of a FogFree™ Strip on mask inner layer.



#### ULTRA™ SENSITIVE FOGFREE™ EARLOOP MASK with Secure Fit™ Mask Technology

CODE	DETAILS	QTY
<b>GCFCXSFSF</b>	White	40/Box

### ULTRA™ SENSITIVE FOGFREE™ EARLOOP MASK WITH SHIELD with Secure Fit™ Mask Technology

- Provides all the features of Ultra™ Sensitive FogFree™ Earloop Mask with Secure Fit™ Mask Technology with addition of the attached optically clear, distortion-free shield.



#### ULTRA™ SENSITIVE FOGFREE™ EARLOOP MASK WITH SHIELD with Secure Fit™ Mask Technology

CODE	DETAILS	QTY
<b>GCPWSSF</b>	White	25/Box

### ULTRA™ EARLOOP MASK with Secure Fit™ Mask Technology

- Fluid resistant outer layer.
- Fluid resistant white inner layer.
- Extra-long aluminum nose piece.



#### ULTRA™ EARLOOP MASK with Secure Fit™ Mask Technology

CODE	DETAILS	QTY
<b>GCFCXUSF</b>	Blue	50/Box

### ULTRA™ SENSITIVE EARLOOP MASK treated with BIOSAFE®

- Fluid resistant outer layer treated with BIOSAFE®, a long-lasting protective bactericidal coating that is permanently bound to the mask and is effective against many harmful bacteria.
- Fluid resistant, extra soft white hypoallergenic inner cellulose layer – will not lint, tear or shred.
- Free inks and dyes.
- Extra-long aluminum nose piece.



#### ULTRA™ SENSITIVE EARLOOP MASK treated with BIOSAFE®

CODE	DETAILS	QTY
<b>GCFCXSBS</b>	White	40/Box

# ULTRA™ SENSITIVE, ULTRA™, BIOSAFE® MASKS CROSSTEX

## TYPE IIR (Enhanced Fluid Resistance)

Ideal for procedures where heavy to moderate amounts of fluid, spray and/or aerosols are produced.

### ULTRA™ SENSITIVE EARLOOP MASK

- Fluid resistant outer layer.
- Fluid resistant, extra soft white hypoallergenic inner cellulose layer – will not lint, tear or shred.
- Free of chemicals, inks and dyes.
- Extra-long aluminum nose piece.



#### ULTRA™ SENSITIVE EARLOOP MASK

CODE	DETAILS	QTY
GCFCXS	White	50/Box

### ULTRA™ SENSITIVE FOGFREE™ EARLOOP MASK

- Provides all the features of Ultra™ Sensitive Earloop Mask with the addition of a FogFree™ Strip on mask inner layer.

#### ULTRA™ SENSITIVE FOGFREE™ EARLOOP MASK

CODE	DETAILS	QTY
GCFCXSF	White	40/Box

### ULTRA™ SENSITIVE FOGFREE™ EARLOOP MASK WITH SHIELD

- Provides all the features of Ultra™ Sensitive FogFree™ Earloop Mask with the addition of attached optically clear, distortion-free face shield.

#### ULTRA™ SENSITIVE FOGFREE™ EARLOOP MASK WITH SHIELD

CODE	DETAILS	QTY
GCPWS	White	25/Box

### ULTRA™ EARLOOP MASK

- Fluid resistant outer layer.
- Fluid resistant white inner layer.
- Extra-long aluminum nose piece.



#### ULTRA™ EARLOOP MASK

CODE	DETAILS	QTY
GCFCXU	Blue	40/Box

### ULTRA™ FOGFREE™ EARLOOP MASK

- Provides all the features of Ultra™ Earloop Mask with the addition of a FogFree™ Strip on mask inner layer.

#### ULTRA™ FOGFREE™ EARLOOP MASK

CODE	DETAILS	QTY
GCFCX	Blue	40/Box

### ULTRA™ FOGFREE™ EARLOOP MASK WITH SHIELD

- Provides all the features of Ultra™ FogFree™ Earloop Mask with the addition of attached optically clear, distortion-free face shield.

#### ULTRA™ FOGFREE™ EARLOOP MASK WITH SHIELD

CODE	DETAILS	QTY
GCPW	Blue	25/Box



## PROCEDURAL MASKS CROSSTEX

### TYPE IIR

Ideal for procedures where moderate to light amounts of fluid, spray and/ or aerosols are produced.

#### PROCEDURAL EARLOOP MASKS

with **Secure Fit™** Mask Technology

- Fluid resistant outer layer.
- Fluid resistant white inner layer.
- Extra-long aluminum nose piece.



**NEW**  
Crosstex Masks  
now also Made  
in Europe

#### PROCEDURAL EARLOOP MASK with **Secure Fit™** Mask Technology

CODE	DETAILS	QTY
<b>GCPBLSFEU</b>	Blue <b>Made in Italy</b>	50/Box
<b>GCPPLVSF</b>	Lavender	50/Box
<b>GCPPKSF</b>	Pink	50/Box

#### PROCEDURAL MEDICAL FACE MASKS

- Fluid resistant outer layer
- Fluid resistant white inner layer
- Extra-long aluminium nose piece
- Comfortable material and ear-loop fit



**NEW**  
Crosstex Masks  
now also Made  
in Europe

#### PROCEDURAL MEDICAL FACE MASKS

CODE	DETAILS	QTY
<b>GCPBLEU</b>	Blue <b>Made in Italy</b>	50/Box
<b>GCPBL100</b>	Blue <b>Made in Italy</b>	100/Box
<b>GCPLV</b>	Lavender	50/Box
<b>GCPPK</b>	Pink	50/Box
<b>GCPYE</b>	Yellow	50/Box

#### SURGICAL TIE-ON MASK

- Fluid resistant outer layer.
- Fluid resistant white inner layer.
- Extra-large panel for surgical procedures.
- Extra-long aluminum nose piece.
- Comfortable material

Made outside of the USA.



#### SURGICAL TIE-ON MASK

CODE	DETAILS	QTY
<b>GCS</b>	Blue	50/Box

# ISOFLUID™ MASKS CROSSTEX

## TYPE II

Ideal for procedures where light amounts of fluid, spray and/or aerosols are produced

### ISOFLUID™ PLUS PURE EARLOOP MASK

with **Secure Fit™** Mask Technology

- Fluid resistant outer layer.
- Fluid resistant, extra soft white hypoallergenic inner cellulose layer – will not lint, tear or shred.
- Free of chemicals, inks and dyes.
- Extra-long aluminum nose piece.



#### ISOFLUID™ PLUS PURE EARLOOP MASK

with **Secure Fit™** Mask Technology

CODE	DETAILS	QTY
<b>GPLUSWHSF</b>	White	50/Box

### ISOFLUID™ EARLOOP MASK

with **Secure Fit™** Mask Technology

- Fluid resistant outer layer.
- White medical grade tissue inner layer.
- Extra-long aluminum nose piece.



#### ISOFLUID™ EARLOOP MASK

with **Secure Fit™** Mask Technology

CODE	DETAILS	QTY
<b>GCIBLSF</b>	Blue	50/Box
<b>GCILVSF</b>	Lavender	50/Box
<b>GCIPKSF</b>	Pink	50/Box

### ISOFLUID™ FOGFREE™ EARLOOP MASK

with **Secure Fit™** Mask Technology

- Provides all the features of Isofluid™ Earloop Mask with Secure Fit™ Mask Technology with the addition of a FogFree™ Strip on mask inner layer.

### ISOFLUID™ FOGFREE™ EARLOOP MASK WITH SHIELD

with **Secure Fit™** Mask Technology

- Provides all features of Isofluid™ FogFree™ Earloop Mask with Secure Fit™ Mask Technology with the addition of attached optically clear, distortion-free face shield.

#### ISOFLUID™ FOGFREE™ EARLOOP MASK

with **Secure Fit™** Mask Technology

CODE	DETAILS	QTY
<b>GCICXBSF</b>	Blue	40/Box

#### ISOFLUID™ FOGFREE™ EARLOOP MASK WITH SHIELD

with **Secure Fit™** Mask Technology

CODE	DETAILS	QTY
<b>GCIPWBSF</b>	Blue	25/Box

### ISOFLUID™ PLUS EARLOOP MASK

- Fluid resistant outer layer.
- Fluid resistant, extra soft white hypoallergenic inner cellulose layer – will not lint, tear or shred.
- Extra-long aluminum nose piece.



#### ISOFLUID™ PLUS PURE EARLOOP MASK

with **Secure Fit™** Mask Technology

CODE	DETAILS	QTY
<b>GPLUSBL</b>	Blue	50/Box
<b>GPLUSKA</b>	Kaleidoscope	50/Box
<b>GPLUSPK</b>	Pink	50/Box

## ISOFLUID™ MASKS CROSSTEX

### TYPE II

Ideal for procedures where light amounts of fluid, spray and/or aerosols are produced

#### ISOFLUID™ EARLOOP MASK

- Fluid resistant outer layer.
- White medical grade tissue inner layer.
- Extra-long aluminum nose piece.



#### ISOFLUID™ FOGFREE™ EARLOOP MASK

- Provides all the features of Isofluid™ Earloop Mask with the addition of a FogFree™ Strip on mask inner layer.

#### ISOFLUID™ FOGFREE™ EARLOOP MASK WITH SHIELD

- Provides all features of Isofluid™ FogFree™ Earloop Mask with the addition of attached optically clear, distortion-free face shield.

##### ISOFLUID™ EARLOOP MASK

CODE	DETAILS	QTY
<b>GCIBL</b>	Blue	50/Box
<b>GCIGR</b>	Green	50/Box
<b>GCILV</b>	Lavender	50/Box
<b>GCIPK</b>	Pink	50/Box
<b>GCISA</b>	Sapphire	50/Box
<b>GCITE</b>	Teal	50/Box
<b>GCITQ</b>	Turquoise	50/Box
<b>GCIBL100</b>	Blue, bulk pack	100/Box

##### ISOFLUID™ FOGFREE™ EARLOOP MASK WITH SHIELD

CODE	DETAILS	QTY
<b>GCIPWB</b>	Blue	25/Box
<b>GCIPWS</b>	Sapphire	25/Box
<b>GCIPWT</b>	Turquoise	25/Box

##### ISOFLUID™ FOGFREE™ EARLOOP MASK

CODE	DETAILS	QTY
<b>GCICXB</b>	Blue	40/Box
<b>GCICXS</b>	Sapphire	40/Box
<b>GCICXT</b>	Turquoise	40/Box

#### LE PETIT MASQUE

- Fluid resistant pink outer layer.
- Soft white spunbound inner layer.
- Smaller than traditional earloop masks, for a better fit.

Made outside of the USA.



##### LE PETIT MASQUE

CODE	DETAILS	QTY
<b>GCAPTPK</b>	Pink	50/Box

## ADDITIONAL FACE PROTECTION CROSSTEX

### ISOLITE™ EARLOOP MASK

- Outer medical grade tissue layer.
- White medical grade tissue inner layer.
- Extra-long aluminum nose piece.

*Does not meet EN14683 nor ASTM classifications.*



#### ISOLITE™ EARLOOP MASK

CODE	DETAILS	QTY
<b>GCLBL</b>	Blue	50/Box
<b>GCLPK</b>	Pink	50/Box

### PATIENT SAFETY MASK WITH SHIELD

- Fluid resistant.
- Covers patient's eyes and nose.
- Protects fragile eye and mucosal tissue from flying particles, spray and splatter.
- Makes patients feel safe, but not "in the dark," with attached clear-vision shield.
- Form-fitting profile, flexible materials provide maximum access to mouth.
- White spunbond inner and outer layers.

*Not made with natural rubber latex.*



#### PATIENT SAFETY MASK WITH SHIELD

CODE	DETAILS	QTY
<b>GCPAT</b>	White	25/Box

### CLEAR FACE SHIELD

- Optically clear, distortion-free wrap-around face shield.
- 1.5 in (3.8 cm) foam headband holds shield away from face; "floats" lightly on forehead, with no pressure on temples; vented for increased air flow.
- Protects mask and face from direct splatter; may prolong mask life.
- Sonically welded elastic headband for added strength.
- Anti-fog treatment on inside and outside of shield.



#### CLEAR FACE SHIELD

CODE	DETAILS	QTY
<b>GCSS</b>	¾ length 17.8 cm	18/Box
<b>GCSS9</b>	Full Length (9") Clear 23 cm	18/Box



## MAXIMUM FILTRATION MASK CROSSTEX

INDICATED  
FOR USE WHEN  
**TREATING  
PATIENTS WITH  
AIRBORNE  
DISEASES.**



### ISOLATOR™ PLUS N95 PARTICULATE RESPIRATOR

- NIOSH approved N95 Respirator.
- Meets guidelines for TB exposure control.
- Resistance to penetration by synthetic blood: 160 mmHg.
- Inner/outer fluid resistant layers.
- Chin pocket construction.
- EN 149: 2001 FFP2 NR / NIOSH N95

*Although the N95 closely resembles a surgical mask, it is a respirator and must be used in accordance with all OSHA regulations regarding respiratory protection.*

CODE	COLOR	PCS
<b>GPRN95</b>	Blue/White Stripe	28/Box

SURGICAL LINE FACE PROTECTIONS AND SHIELDS OMNIA



30.Z1025.00

**Omniprotector** is a very light protective shield mounted on a PVC frame. Its optical Polyethylene lenses offer a particularly clear viewing and do not alter images. Available in 2 models, one with nose rest.

CODE	DESCRIPTION	QTY
30.Z1025.00	Protection eyewear with disposable visor: the package includes 3 supports and 20 visors	1/Box
30.Z1045.00	Protective glass with nose clip and disposable visor, 1 support and 20 disposable visors	1/Box
30.Z1026.00	Spare visor for the 30.Z1025 and 30.Z1045 models	20/Box




30.Z1045.00



10.M0200.00

The **Visual** is the best individual protection device available to ensure your safety. Its wide polyethylene fog preventing shield and the polyurethane insert make this device also fit for use by people wearing glasses.

CODE	QTY
10.M0200.00	12/Box



30.Z1040.00

Protection eyewear with nose clip and long disposable visor. The package includes 1 PVC support and 20 spare visors.

CODE	QTY
30.Z1040.00	1/Box

Spare visor for the 30.Z1040 model.

CODE	QTY
30.Z1039.00	20/Box



30.M0245.00

**Omnishield** is a light, transparent plastic eyewear with side shield designed to ensure excellent protection without interfering with the visual field. This accessory is absolutely necessary in dental procedures.

CODE	QTY
30.M0245.00	5/Box

For more information on Hu-Friedy products or to place an order, call at 00800 48 37 43 39 or visit [HuFriedyGroup.eu](https://HuFriedyGroup.eu)

For more information on Crosstex products, call +39 05 24 52 74 53 or visit [Crosstex.com](https://Crosstex.com)

For more information on Omnia products or to place an order, call +39 05 24 52 74 53 or visit [OmniaSpa.eu](https://OmniaSpa.eu)

Hu-Friedy Mfg. Co., LLC • European Headquarters • Lyoner Str. 9 • 60528 Frankfurt am Main, Germany • [HuFriedyGroup.eu](https://HuFriedyGroup.eu)

As some of the products included in this catalogue might not be available in some countries, we invite you to verify the availability of the product with your local authorised dealer or to contact our customer care at the free phone number 00800 48 37 43 39. Available while stocks last. All other company and product names are trademarks of Hu-Friedy Mfg. Co. LLC, its affiliates or related companies, unless otherwise noted.

©2021 Hu-Friedy Mfg. Co., LLC All rights reserved. 90.L0390.EN/1221