



HU-FRIEDY SIGNATURE SERIES PROCEDURE TUBS

GET ORGANIZED & BE EFFICIENT



SIGNATURE SERIES **PROCEDURE TUB BENEFITS:**

STANDARDIZES MATERIAL SET-UP FOR EVERY TYPE OF PROCEDURE

- Increases efficiency of materials management
- Eliminates time-consuming tray preparation for every procedure
- Tubs may be configured to meet your office needs or are available as a complete tub

PROCEDURE TUB COLORS MATCH HU-FRIEDY IMS CASSETTE RAIL COLORS

- Increases office organization and workflow efficiency
- Allows for guick identification of consumable products
- Provides esthetic integration with IMS cassettes

SIGNATURE SERIES TUBS

LOCKING COVERS

- Protects items from airborne materials during transportation and storage
- Adheres to infection control standards

ANTIMICROBIAL PROTECTION

- Inhibits the growth of microorganisms that can cause stains, odors and product degradation
- Registered with the EPA for these applications

ALLOWS FOR STORAGE IN CENTRAL STERILIZATION OR EACH OPERATORY

• Tubs can be transported to the operatory prior to each procedure, eliminating inventory management of materials in each operatory

ANTIMICROBIAL PROTECTION

Signature Series Tub

Tub Divider (package of 8)

• Single Tub Cup w/Cover*

Double Tub Cup w/Cover*

Divided Slide Tray

Locking Tub Cover*

Tub Liner*

Antimicrobial protection is integrated into the tub and tub components during the manufacturing process. This technology works by interrupting key cell functions that allow the microorganisms to function, grow and reproduce. Antimicrobial protection begins to work as soon as the microorganism comes into contact with the product surface and works continuously to inhibit the growth of microbes that can cause stain, odors and product degradation - for the lifetime of the product.

22 4 y 20 4 y 10 2 cm

SIGNATURE SERIES COMPLETE TUB

	32,4 x 28,4 x 10,2 cm	
	Red	IMS-1400
	Gray	IMS-1401
	White	IMS-1402
	Orange	IMS-1403
	Purple	IMS-1404
	Yellow	IMS-1405
	Lavender	IMS-1406
	Black	IMS-1407
	Blue	IMS-1408
	Green	IMS-1409
	Ocean Blue	IMS-14011

antimicrobial protection.

SIGNATURE SERIES TUB ACCESSORIES

SIGNATURE SERIES TUB

1 x 27,6 x7cm	
ed	IMS-1410
ay	IMS-1411
hite	IMS-1412
ange	IMS-1413
ırple	IMS-1414
llow	IMS-1415
vender	IMS-1416
ack	IMS-1417
ue	IMS-1418
een	IMS-1419
cean Blue	IMS-14111



DIVIDED SLIDE TRAY*

24,8 x 13,7 x 2,7cm IMS-1421

 Can be used to store cotton rolls, cotton squares, clamps, syringe tips, etc.



FLAT SLIDE TRAY*

24,8 x 13,7 x 2,7cm IMS-1422

 Can be used to store discs. strips, toffelmeil, articulating paper, etc.



LOCKING TUB COVER*

32,4 x 28,4 x 7cm IMS-1420

 Protects items from airborne materials during transportation and storage



TUB INSERT

24,8 x 27,3 x 2,6cm IMS-1423

- The insert features 10 compartments for easy organization of tub materials
- Sits on rim of tub allowing for storage underneath



Includes 8 dividers

TUB LINER*

DIVIDERS

24,5 x 22,5 x 0,3cm IMS-1428

Non-slip liner prevents small items from slipping under divider sections

IMS-1429

IMS-1430

SIGNATURE SERIES CLIP ON HOLDERS





CAPSULE CLIP-ON HOLDER

20,3 x 6,4 x 1,9cm

• Allows quick access to composite material during procedure

Conveniently clips on any side of tub

• Holds up to 8 capsules, 4 syringe tubes and capsule/gun dispenser

SYRINGE CLIP-ON HOLDER

20,3 x 6,7 x 1,9cm

- Allows quick access to composite material during procedure Conveniently clips on any side of tub
- Holds up to 15 syringe tubes and 3 applicator tips

SIGNATURE SERIES TUB CUPS AND COVERS

Can be used to store brushes, clamps, cotton pellets, syringe tips, prophy SINGLE TUB WITH COVER* paste, prophy angle tips, pumice, implant tools, etc.



5,7 x 4,1 x 3cm	IMS-1424

DOUBLE TUB WITH COVER*

8,6 x 5,7 x 3cm	IMS-1425
-----------------	----------

LONG TUB WITH COVER*

21,2 x 5,7 x 3cm

IMS-1426

- Provides 5 storage compartments
- Dividers are removable for customization
- · Perfect for storing capsules, discs, abrasive strips, wedges, matrix bands, etc.

SIGNATURE SERIES SYRINGE STANDS



20-UNIT SYRINGE STAND

22.9	x 12	1 x	6.4	lcm	

- Stands up for easy access to syringes
- · Lays flat to fit inside a tub or drawer
- Holds up to 20 syringe tubes

10-UNIT SYRINGE ADD ON

21,6 x 6 x 1,9cm

IMS-1432

IMS-1431

- An add-on for the Syringe 20 Unit Stand
- Provides additional space for syringe storage
- Holds up to 10 syringe tubes

SIGNATURE SERIES RACKS





Adiustable

ADJUSTABLE 3 TUB RACK

MIN 19,05 cm - MAX 39,37 x 33,02 x 26,67 cm IMS-1433

- · Holds a maximum of 3 tubs with lids or 6 cassettes
- Non-slip feet
- Sturdy metal construction
- Includes fasteners to secure and stabilize
- No screws, tools, nuts or bolts needed

ADJUSTABLE 4 TUB RACK

MIN 19,05 cm - MAX 39,37 x 26,67 x 47,62 cm IMS-1434

- · Holds a maximum of 4 tubs with lids or 8 cassettes Non-slip feet
- Sturdy metal construction
- Includes fasteners to secure and stabilize
- No screws, tools, nuts or bolts needed

IMS CASSETTE AND TUB SYSTEM

HOW CASSETTES AND TUBS WORK TOGETHER:

- Together, procedure tubs and cassettes create the most efficient and organized system for managing instruments and consumable products.
- Time savings achieved by using IMS (Instrument Management System) cassettes and procedure tubs will allow the staff to spend time on more value added activities in the dental office.
- Procedure tubs can be used as portable operatory drawers to store, organize and transport consumable materials to and from the storage/sterilization area and the operatory.
- All of the materials needed for a specific procedure can be stored in the tub. A different tub and color can be used for each different procedure type (crown & bridge, composite, etc.).
- The Signature Series procedure tubs are available in colors that match IMS cassette rails, which allow the staff to match the tubs and cassettes by procedure.

Sealants

PROCEDURE TUB SET-UP

EXAMPLES



Crown and Bridge

AN EXAMPLE OF HOW IT WORKS:

Endontics

Composite

• Blue cassette rails and a blue tub may represent a composite procedure. All of the consumable materials for composite procedures can be stored and organized in the blue tub and instrumentation stored in the cassette with blue rails, allowing for quick identification.

ORDERING INFORMATION

Product Description Part Codes IMS-1400 Complete Tub - Red Complete Tub - Gray IMS-1401 Complete Tub - White IMS-1402 Complete Tub - Orange IMS-1403 Complete Tub - Purple IMS-1404 Complete Tub - Yellow IMS-1405 IMS-1406 Complete Tub - Lavender Complete Tub - Black IMS-1407 Complete Tub - Blue IMS-1408 Complete Tub - Green IMS-1409 Complete Tub - Ocean Blue IMS-14011 Tub Only - Red IMS-1410 Tub Only - Gray IMS-1411 Tub Only - White IMS-1412 IMS-1413 Tub Only - Orange IMS-1414 Tub Only - Purple IMS-1415 Tub Only - Yellow IMS-1416 Tub Only - Lavender

	Product Description	Part Codes
	Tub Only - Black	IMS-1417
	Tub Only - Blue	IMS-1418
	Tub Only - Green	IMS-1419
_	Tub Only - Ocean Blue	IMS-14111
	Locking Tub Cover	IMS-1420
	Divided SlideTray	IMS-1421
	Flat Slide Tray	IMS-1422
_	Tub Insert	IMS-1423
	Tub Cup Single With Cover	IMS-1424
_	Tub Cup Double With Cover	IMS-1425
	Tub Cup Long With Cover	IMS-1426
_	Tub Dividers	IMS-1427
	Tub Liner	IMS-1428
_	Capsule Clip-On Holder	IMS-1429
	Syringe Clip-On Holder	IMS-1430
_	20-Unit Syringe Stand	IMS-1431
	10-Unit Syringe Add On	IMS-1432
_	Adjustable 3 Tub Rack	IMS-1433
	Adjustable 4 Tub Rack	IMS-1434

(Surface disinfect only. Do not heat sterilize Signature Series Tubs and Accessories.)

For more information about Hu-Friedy's Total Solutions visit

HU-FRIEDY.EU









INSTRUMENT

MANAGEMENT







ULTRASONIC

SCALING

Manufacturer: Hu-Friedy Mfg. Co., LLC | 3232 N. Rockwell Street | Chicago, IL 60618 | USA Customer Care Department:Free Phone 00800 48 37 43 39 | Free Fax 00800 48 37 43 40 Website: Hu-Friedy.eu



©2014 Hu-Friedy Mfg. Co., LLC. All rights reserved, IMS-380GB/0914