


Diamond Dusted Forceps

Get a Better Grip



Our **Diamond Dusted Forceps** feature a fine diamond coating which increases resistance between the forceps and tooth surface providing enhanced engagement. The secure grasp of the crown and root reduces the risk of crushing or fracturing the tooth.

A light diamond coating bonded to longitudinal beak serration to significantly increase the grip compared to traditional forceps.*



Thin, tapered beak designs allow for improved access compared to traditional forceps beaks.







Matte finish for additional contrast and reduced light reflection.




Unique hole design reduces the weight of the forceps, providing increased comfort and control.

Diamond Dusted Forceps





are available in the following Atraumair patterns:




Upper

F1XSD	1 Standard Forceps, Diamond Dusted Upper Incisors	
F35XSD	35 Forceps, Diamond Dusted Upper Canines and Premolars	
FX44XSD	44 Root Forceps, Diamond Dusted Upper Roots	
F67XSD	67 Apical Forceps, Diamond Dusted Upper Molars	

FAF1XSD	1 Apical Forceps, Diamond Dusted Upper Incisors	
FAF35XSD	35 Apical Forceps, Diamond Dusted Upper Canines and Premolars	
FX51XSD	51 Root Forceps, Diamond Dusted Upper Roots	

Lower




F13XSD	13 Forceps, Diamond Dusted Lower Premolars	
F36XSD	36 Forceps, Diamond Dusted Lower Incisors	
FX45XSD	45 Root Forceps, Diamond Dusted Lower Roots, Serrated	
FX79XSD	79 Forceps, Diamond Dusted Lower Molars	

F21XSD	21 Forceps, Diamond Dusted 1st and 2nd Molars	
FAF36XSD	36 Apical Forceps, Diamond Dusted Lower Premolars and Incisors	
FX73XSD	73 Forceps, Diamond Dusted Lower Molars	

Find your optimal extraction forceps at **HuFriedyGroup.eu**

Free Call: 00800 48 37 43 39 | Free Fax: 00800 48 37 43 40

E-Mail: info@hu-friedy.eu | Website: HuFriedyGroup.eu

   HuFriedyGroupEurope

* Increased grip by 20% compared to standard Hu-Friedy Forceps. Data on file.
Hu-Friedy Mfg. Co., LLC • European Headquarters • Lyoner Str. 9 • 60528 Frankfurt am Main, Germany • HuFriedyGroup.eu
All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted.
©2023 Hu-Friedy Mfg. Co., LLC. All rights reserved. HFL-666GB/0223