

# HuFriedyGroup

## Key Opinion Leader

### Dr. Alfonso Rao

Instrument list  
& Omnia Set for Implantology



Focus on:  
**Implant Surgery**



Country: **UK**



Favourite Instrument:  
**Buser Periosteal Elevator  
PPBUSER6**



Website:  
**[www.deltadentalacademy.co.uk](http://www.deltadentalacademy.co.uk)**



Inspired by:  
**Istvan Urban,  
Egon Ewue**




Quality instruments are essential to achieve the best outcome, with HuFriedyGroup the quality is assured!

Dr. Alfonso Rao




## DR. ALFONSO RAO INSTRUMENT LIST

This kit enables a broad range of treatment procedures to be performed using the minimum number of ergonomically designed dental implant surgery instruments. Each instrument is engineered to guarantee ease of use. The quality of your instrument can vastly improve your experience and for this reason I have choose HuFriedyGroup instruments for my practice and for Delta Dental Academy. HuFriedyGroup instruments unique quality, efficiency and reliability provides you with the most advanced dental product solutions when learning your craft and improving your skills.



Visit us online at **[HuFriedyGroup.eu](http://HuFriedyGroup.eu)**

# INSTRUMENT LIST




**10-130-05E**  
Scalpel Handle #5E,  
European handle


**PPBUSER6**  
Periosteal Buser hdl #6

**PJO/16**  
Chisel Paro Jovanovic  
#1 hdl #6 4 mm

**CL866**  
Surgical Curette Lucas  
#86, hdl #6 2,8 mm



**H4**  
Hemostat Halsted-Mosquito  
#4 straight, 12 cm




**NH5052**  
Fine Swedish Pattern Perma Sharp™  
Needle Holder, 5-0, 6-0, 7-0, 8-0  
15 cm

**PERMA SHARP™  
INSTRUMENTS**  
Featuring tungsten carbide  
inserts bonded to surgical-  
grade stainless steel on both  
cutting edges for smooth,  
controlled cutting and  
maximum durability.



**MHE6**  
Mirror Handle #6  
European Simple Stem





**MIR5HDE**  
#5 HD Single-Sided Mirror,  
EU Simple Stem




**CRM**  
Cheek Retractor Minnesota

**SUPER-CUT**  
Cutting blade honed/  
serrated to „razor“  
sharpness with opposing  
blade finished with fine  
serrations to avoid  
tissue slippage.



**SCORTSC**  
Super-Cut Scissors  
Cortellini curved, 12 cm



**TPG3**  
Gerald 1x2 Tissue Pliers, Straight, 18 cm




**TP5050**  
Perma Sharp™ Tissue Pliers Gerald, straight point, 18 cm



## INSTRUMENT LIST



**PPR36**  
Periosteal Prichard #3 hdl #6



**PCPUNC156**  
Probe UNC #UNC15 hdl #6



**CSEUWES6**  
Euwe Bone Scoop/Plugger hdl #6

### CHOOSE THE COLOR.

Just add the color coding number to the last part of the partcode (replace " " with the number)

Red	0
Grey	1
White	2
Orange	3
Purple	4
Yellow	5
Lavender	6
Black	7
Blue	8
Green	9
Ocean Blue	11

### WHY YOU SHOULD USE AN INSTRUMENT MANAGEMENT SYSTEM (IMS):

By standardising your instruments with IMS System, additional time and revenue is generated while extending the useful life of your Hu-Friedy instruments. Dental professionals state that by implementing IMS™ Cassettes, they save 5 to 10 minutes per procedure allowing the staff to spend more time on value-added activities\*. In addition, safer transportation of soiled instruments reduces risk of injury and exposure to blood borne pathogens.

\*When compared to single instrument decontamination. Based on market survey results.

DISCOVER MORE AT IMS WEBPAGE



### IMEDIN18\_S

IMS™ Infinity Series™ **DIN Surgical Cassette for 18 instruments, with an accessory area + 4 clips (2 x I1000) for hinged instruments.**  
Dimensions: 290 mm x 184 mm x 34 mm



## ALFONSO RAO'S OMNIA SET FOR IMPLANTOLOGY



This Kit has been designed for experienced clinicians and for those who want to perform complex surgery especially in the Aesthetic zone. The Microsurgical and minimally invasive surgery make it possible to minimise the tissue trauma and significantly enhance the final outcome.

Dr. Alfonso Rao



### 12.E1013.00

Alfonso Rao's Omnia Set for Implantology

Pack contains irrigation lines compatible with most models of implantology motors, for other further information please contact your HuFriedyGroup representative.

### DESCRIPTION

- 2 x Special Plus model gown (length 115 cm) size L
- 2 x Round cap with elastic
- 1 x Surgical cap with fixed ties
- 2 x 3-ply mask with round anallergic elastics light blue
- 1 x Light blue absorbent/waterproof drape 50 x 50 cm
- 2 x Light blue absorbent/waterproof drape 50 x 75 cm
- 1 x Absorbent/Waterproof drape 50 x 75 cm with adhesive U-shaped split 11 x 9 cm
- 5 x Gauze 10 x 10 cm with 4 layers, White NWF
- 1 x Sheaths 120 x 7 cm with inserter Omnisleeve
- 5 x Adhesive film sheets 20 x 20 cm 2Cross20
- 1 x Irrigation set - (partcode 32.F0134.00 - compatible with the following micromotors: NSK Surgic XT, Surgic XT Plus, Osseocision®, VarioSurg, SurgicPro, Kyocera Implantor, Neo, Bonart Piezo surgical system, Dentium iCTmotor, NSK 3i Osseocision, Delfine Intra-Lock)
- 1 x Cannula, model OmniaSurg
- 1 x Transparent bag 67 x 75 cm with sliding ties for physiodispensers
- 1 x Bowl 250 ml PE
- 1 x Dish 250 ml + patient disinfection swab
- 1 x Bag with adhesive tape for collecting contaminated waste



Schedule a virtual consultation to learn how IMS Cycle can prevent the spread of disease and improve the safety of your practice.



**Chris Mason (Country Manager UK & Ireland)**  
Tel.: 07703 186 474 | cmason@hu-friedy.com



## DELTA DENTAL ACADEMY



### Practice-ready dental concepts from world-leading trainers

Delta Dental Academy is a dynamic multidisciplinary training centre in the heart of Bristol. Delegates learn advanced techniques from some of the world's top dentists sharing a groundbreaking syllabus that's designed to be ready to put into practice from the moment you leave the classroom. Delta Dental Academy is one of the UK's few training centres to offer expert-led training on patients provided by the academy. Our unique blend of theory and practical training has empowered hundreds of dentists to learn, grow and develop their practice.



#### Get Practice-Ready Knowledge

Our syllabus is designed for dentists to put theory into practice from the moment they leave the classroom.



#### Learn with Real Patients

Delta Dental Academy is one of the UK's few dental academies to provide delegates with live patients to train on.



#### Expert-Led Training

Our faculty comprises some of the leading experts and educators in endodontics, implantology, oral surgery and more.



#### Top-Rated by Dentists

With hundreds of happy cohorts trained over almost two decades at the forefront of dental education, Delta is consistently top-rated by dental professionals.

## DELTA DENTAL ACADEMY

### COURSES

- 10-Day Endodontics
- The Famous 4 on the case
- Occlusion and Restorative Dentistry for the GDP
- Nurses introduction to dental Implants
- Introduction to Soft Tissue Management
- Guided Bone Regeneration and Placement in the Aesthetic Zone
- Endodontic Retreatment Masterclass
- Full-Arch: Immediate Placement & Loading
- Implant Fundamentals
- Essential Endodontics - A Modern Approach for GDPs
- Micro-surgical Endodontics
- Master in Microendodontics

### ONLINE COURSES

- Endodontic & Restorative Management of Complex Cases
- Micro Apical Placing of MTA into Wides Apex
- Master MTA




[www.deltadentalacademy.co.uk](http://www.deltadentalacademy.co.uk)



**For more information or to place an order, call at 00800 48 37 43 39 or visit [HuFriedyGroup.eu](https://www.HuFriedyGroup.eu)**

**Free Call:** 00800 48 37 43 39 | **Free Fax:** 00800 48 37 43 40

**E-Mail:** [info@hu-friedy.eu](mailto:info@hu-friedy.eu) | **Website:** [HuFriedyGroup.eu](https://www.HuFriedyGroup.eu)

   [HuFriedyGroupEurope](https://www.HuFriedyGroupEurope.com)

Hu-Friedy Mfg. Co., LLC • European Headquarters • Lyoner Str. 9 • 60528 Frankfurt am Main, Germany • [HuFriedyGroup.eu](https://www.HuFriedyGroup.eu)  
All company and product names are trademarks of Hu-Friedy Mfg. Co. LLC., its affiliates or related companies, unless otherwise noted.  
©2023 Hu-Friedy Mfg. Co., LLC. All rights reserved. HFL-616GB/0323